

# Can a patient wear their glasses or contact lenses during visual field testing on a Matrix 800?

## ZEISS Quick Help: Humphrey® Matrix 800

Yes, patients can wear their normal glasses, contacts or a trial frame if needed during visual field testing on a Matrix 800. It is okay to use bi-focal or progressive lenses. Use a trial frame instead of the patient's normal glasses or contacts if the patient's corrective lenses are tinted or photo-chromatic, or if the patient's frame obscures part of the testing field.


For accurate visual field test results, correct the patient to within the refractive error cited in the table below. Verify that the patient's eyes are not dilated during the visual field test, unless the pupil is less than 3 mm.

Tests patient corrected within:	Tests patient corrected within:
N-30-5 (-1) FDT Screening N-30-F FDT Threshold	+/- 6 Diopters
24-2-5 (-1) FDT Screening 24-2 FDT Threshold 30-2 FDT Threshold	+/- 3 Diopters
Macula FDT Threshold 10-2 FDT Threshold	+/- 2 Diopters

Refer to the Humphrey Matrix user manual *Instructions for Use* for safe and effective operation of the instrument.



 **Carl Zeiss Meditec AG**  
Goeschwitzer Strasse 51–52  
07745 Jena  
Germany  
[www.zeiss.com/med](http://www.zeiss.com/med)

 **Carl Zeiss Meditec, Inc.**  
5160 Hacienda Drive  
Dublin, CA 94568  
USA  
[www.zeiss.com/med](http://www.zeiss.com/med)



**SER.9255** Printed in the United States. CZ-1W2017 United States edition: Only for sales in selected countries. The contents of the reference guide may differ from the current status of approval of the product or service offering in your country. Please contact our regional representatives for more information. Subject to change in design and scope of delivery and as a result of ongoing technical development. Humphrey is either a trademark or registered trademark Carl Zeiss Meditec AG or other companies of the ZEISS Group in Germany and/or other countries. © Carl Zeiss Meditec, Inc., 2017. All rights reserved.