


## Changing the Zeiss user password for the NAS device




### Note

When changing the Zeiss user password for the NAS, it must match the Zeiss user password for the Windows login.

- 1) Verify the NAS device is powered on.
- 2) Log out and exit the PLEX Elite application.
- 3) Launch Internet Explorer.
- 4) In the address field, enter **//kilinas** and then press **Enter** on the keyboard.
- 5) When the NAS login window appears, enter **zeiss** for the user name.
- 6) For the password, enter the initial password that is used by the Zeiss user when logging into the PLEX Elite instrument. (For that initial password, refer to section 6 of the PLEX Elite version 2.0 IFU or contact your ARI ambassador.)
- 7) Select **Sign In**. This will open the Synology DiskStation Manager (DSM).



8) Select the Options icon from the taskbar.



9) Enter the current password, enter the new password and then enter the new password again to confirm.



**Note**

The new password for the NAS must match the Zeiss user password for Windows.

A screenshot of the Windows 'Personal' account settings window. The window has a title bar with 'Personal' and standard window controls. Below the title bar are tabs: 'Account', 'Account Protection', 'Quota', 'Desktop', 'Email Account', and 'Others'. The 'Account' tab is selected. The form contains the following fields:


- Name: zeiss
- Description: (empty text box)
- Password: (empty text box)
- New Password: (text box with 7 dots)
- Confirm password: (text box with 7 dots)
- Email: user@email
- Display language: System setting (dropdown menu)
- Date format: YYYY-MM-dd (System setting) (dropdown menu)
- Time format: 24-hour time (System setting) (dropdown menu)

Below these fields is a checkbox for 'Enable 2-step verification' which is unchecked. There is a '2-Step Verification' button. Below that is a link for 'View your account activity, including current connections, trusted devices, and login history.' and an 'Account Activity' button. At the bottom right are 'OK' and 'Cancel' buttons.

10) Select **OK**.


11) If you wish to enable “Auto Block,” perform the following steps.

- a) Select **Control Panel > Security > Account**.
- b) Check box for “Enable auto block.”
- c) If desired, change any of the Auto Block default settings.
- d) Select **Apply**.



12) If you wish to disable the “Tech Support” user, perform the following steps.

- a) Select **Control Panel > User (or Control Panel > User & Group)**.
- b) Select the Tech Support user.
- c) Select **Edit**.
- d) Check box for “Disable this account.”
- e) Select **Save**.



Disable this account

Immediately

13) Close the DSM window.