Microsoft Windows Print Spooler Vulnerability

CVE-2021-34527

("PrintNightmare") Cybersecurity Update For Windows 7 OS


Table of Contents

English	3
العربية.....	8
български.....	12
中文 (simplified).....	16
Hrvatski	20
Česky	24
Dansk	28
Dutch.....	32
Eesti	36
Suomi	40
Français.....	44
Deutsch	48
ελληνικά.....	52
Magyarul	56
Italiano	60
日本語.....	64
한국어.....	68
Latviešu.....	72
Lietuvių	76
Norsk.....	80
Polski.....	84
Português (Portugal).....	88
Português (Brasil)	92
Român.....	96
Русский	100

Table of Contents

Slovensky.....	104
Slovenščina	108
Español.....	112
Español (Latin America)	116
Svenska.....	120
Türkmen.....	124

Microsoft Windows Print Spooler Vulnerability

CVE-2021-34527

("PrintNightmare") Cybersecurity Update For Windows 7 OS

Copyright

© 2021, Carl Zeiss Meditec, Inc., Dublin, CA

Trademarks

All Zeiss products mentioned herein are either registered trademarks or trademarks of Carl Zeiss Meditec, Inc. in the United States and/or other countries.

All other trademarks used in this document are the property of their respective owners.

Patents

www.zeiss.com/meditec/us/imprint/patents.html

1 About the Update

PrintNightmare is a vulnerability affecting Microsoft Windows operating systems (OS).

A remote code execution vulnerability exists when the Windows Print Spooler service improperly performs privileged file operations. An attacker who successfully exploited this vulnerability could run arbitrary code with SYSTEM privileges. An attacker could then install programs; view, change, or delete data; or create new accounts with full user rights.^[1]

PrintNightmare does not affect safety and performance on any of the ZEISS devices.

However, ZEISS recommends updating devices with the Microsoft patch and/or registry settings, as applicable, to ensure continued cybersecurity. ZEISS has analyzed the impact of the vulnerability on ZEISS products running Windows OS and only devices listed below must run the update.

- CIRRUS 400/4000
- CIRRUS 500/5000
- CIRRUS Photo 600/800
- ATLAS 9000
- HFA3
- PLEX Elite 9000

2 Disable the Printer Spooler Group Policy

It is recommended that a site IT Administrator performs this task.

1. Login to the instrument.
2. In the search box on the taskbar, type `gpedit.msc` to run the Group Policy editor.
3. Navigate to **Local Computer Policy > Computer Configuration > Administrative Templates > Printers**.
4. Select **Allow Print Spooler** to accept client connections.
5. Double-click the policy to open it.
6. Select **Disabled**.
7. Select **OK** to save the policy.
8. Restart the instrument.

Result

- ✓ Doing the above steps will block the remote attack vector by preventing inbound remote printing operations. The system will no longer function as a print server, but local printing to a directly attached device will still be possible.

^[1] Microsoft website: <https://msrc.microsoft.com/update-guide/vulnerability/CVE-2021-34527>

- ✓ This Group Policy change will set the following registry key
`HKEY_LOCAL_MACHINE\ Software\Policies
\Microsoft\Windows NT\Printers\ Register-
SpoolerRemoteRpcEndPoint = 2`

3 Deactivate Point and Print in the Windows Registry

We recommend an IT administrator to perform this task.

Action

1. Login to the instrument.
2. In the search box on the taskbar, type `regedit`, then select Registry Editor.
3. Navigate to `HKEY_LOCAL_MACHINE\SOFTWARE\Policies\Microsoft\Windows NT\Printers`.
4. If it exists, expand the Printers branch, and make sure that the PointandPrint group does not exist.

NOTE! The Printers branch does not exist in the factory configuration.

- ⇒ If the **PointandPrint** group exists, check the following settings, if it exists then set the value to 0 .
- `NoWarningNoElevationOnInstall = 0` or does not exist
 - `UpdatePromptSettings = 0` or does not exist

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