

SL Workstation cybersecurity update ("PrintNightmare")



Seeing beyond

SL Workstation cybersecurity update ("PrintNightmare")

Please note: This document does not replace the user manual which is delivered with the device.

About the Update

"PrintNightmare" is a vulnerability affecting Microsoft Windows operating systems (OS). According to Microsoft, "a remote code execution vulnerability exists when the Windows Print Spooler service improperly performs privileged file operations. An attacker who successfully exploited this vulnerability could run arbitrary code with SYSTEM privileges. An attacker could then install programs; view, change, or delete data; or create new accounts with full user rights."

"PrintNightmare" does not affect the safety and performance of the ZEISS SL Workstation. Nevertheless, ZEISS would like to offer you an update for your ZEISS SL Workstation to fix the "PrintNightmare" vulnerability.


Required storage media

You will need an empty USB flash drive with at least 1 GB of storage capacity to download the update files to prior to installation.


Preparation

- Download the following files to the root directory of your USB flash drive:
 - LGPO.exe
 - PrintNightmare_1.lgpo.txt
 - Install_Patch_PrinterNightmare_CVE-2021-34527_Win-7.cmd


How to update your SL Workstation




1. Power on the SL Workstation.
2. Connect the USB flash drive containing the update files.
3. Copy all three files to the following location
c:\Users\zeiss\Desktop\scripts. Check that the .cmd file name contains "Win-7".




4. Click the start menu icon in the system taskbar and enter "cmd" in the search field. A view similar to the image on the left appears.



5. Right click on "cmd.exe" to open the context menu.
Select "Run as administrator".



6. A new window appears on the PC screen. Select "Yes" to start the command prompt.



7. A command prompt with administrator rights opens.
Enter "cd c:\Users\zeiss\Desktop\Scripts" and press Enter.

```
Administrator: C:\Windows\System32\cmd.exe
Microsoft Windows [Version 6.1.7601]
Copyright (c) 2010 Microsoft Corporation. All rights reserved.

C:\Windows\System32>cd C:\Users\zeiss\Desktop\Scripts
C:\Users\zeiss\Desktop\Scripts>Install_Patch_PrinterNightmare_CVE-2021-34527_Win-7.cmd
```

8. After the directory has been changed, enter "Install" and press the TAB key. The filename should be completed as shown in the image on the left. Press Enter to continue.

```
Administrator: C:\Windows\System32\cmd.exe - Install_Patch_PrinterNightmare_CVE-2021-34527_Win-7.cmd
C:\Users\zeiss\Desktop\Scripts>Install_Patch_PrinterNightmare_CVE-2021-34527_Win-7.cmd
OK : 1 LGPO compiled and registered.

SERVICE_NAME: lpdsvc
        TYPE               : 20  WIN32_SHARE_PROCESS
        STATE                : 3  STOP_PENDING
                        (NOT_STOPPABLE, PAUSABLE, ACCEPTS_SHUTDOWN)
        WIN32_EXIT_CODE       : 0  (0x0)
        SERVICE_EXIT_CODE   : 0  (0x0)
        CHECKPOINT           : 0x0
        WAIT_HINT            : 0x0

SERVICE_NAME: spooler
        TYPE               : 110 WIN32_OWN_PROCESS (interactive)
        STATE                : 3  STOP_PENDING
                        (NOT_STOPPABLE, NOT_PAUSABLE, IGNORES_SHUTDOWN)
        WIN32_EXIT_CODE       : 0  (0x0)
        SERVICE_EXIT_CODE   : 0  (0x0)
        CHECKPOINT           : 0x3
        WAIT_HINT            : 0x4e20
This shared resource does not exist.
More help is available by typing NET HELPMSG 2310.


SERVICE_NAME: spooler
        TYPE               : 110 WIN32_OWN_PROCESS (interactive)
        STATE                : 2  START_PENDING
                        (NOT_STOPPABLE, NOT_PAUSABLE, IGNORES_SHUTDOWN)
        WIN32_EXIT_CODE       : 0  (0x0)
        SERVICE_EXIT_CODE   : 0  (0x0)
        CHECKPOINT           : 0x0
        WAIT_HINT            : 0x7d0
        PID                 : 3588
        FLAGS                :

SERVICE_NAME: lpdsvc
        TYPE               : 20  WIN32_SHARE_PROCESS
        STATE                : 2  START_PENDING
                        (NOT_STOPPABLE, NOT_PAUSABLE, IGNORES_SHUTDOWN)
        WIN32_EXIT_CODE       : 0  (0x0)
        SERVICE_EXIT_CODE   : 0  (0x0)
        CHECKPOINT           : 0x0
        WAIT_HINT            : 0x7d0
        PID                 : 1392
        FLAGS                :


OK : Print spooler service restarted.
Press any key to continue . . .
```

9. The script is executed.


10. After the script has run successfully and completely, restart the computer for the changes to become effective.



11. After the restart, click the start menu icon in the taskbar. Enter "gpedit.msc" in the search field. Click on the item that appears in the list to open the group policy settings.



12. The following message will appear two times. Confirm the message with “OK”.



13. The group policy settings window opens. Go to “Computer Configuration” > “Administrative Templates” > “Printers”.

Check that the setting for “Point and Print Restrictions” is set to “Enabled” and “Allow Print Spooler to accept client connections” is set to “Disabled”.

Discover more expert videos, supporting documents, and common questions and answers in the ZEISS Product Insights.



ZEISS Product Insights website



Carl Zeiss Meditec AG

Goeschwitzer Strasse 51-52

07745 Jena, Germany

www.zeiss.com/med

info.meditec@zeiss.com

000000-2279-200-AddGA-GB-260821