

- Remove the defective fuse (**2**, Fig. 25) from the fuse carrier (**1**, Fig. 25) and replace it with a new fuse.
- Replace the fuse compartment (**1**, Fig. 25) and ensure the locking mechanism clicks in.

## Maintenance

### Care and cleaning



#### **WARNING - RISK OF ELECTRIC SHOCK**

Repair and maintenance work may only be carried out by authorized service staff.

Do not remove any parts of the casing or carry out any repairs yourself, especially not on the interior parts of the device.

To avoid overheating or the risk of fire, only replace spent fuses with new ones of the same type.

Removing parts of the casing exposes you to the risk of an electric shock.

- There are no parts inside the device which can be serviced by the user.
- Any service work which needs to be performed must be carried out by trained personnel.

#### **CAUTION - PROPERTY DAMAGE**

The national disinfecting regulations must be observed in the choice of disinfectants and disinfection procedures. Please note that some cleaning agents and disinfectants may have an adverse effect on plastic components. Damage caused by such disinfectants is not covered by our warranty. The surfaces of the device have been tested and are guaranteed to resist frequent treatment with alcoholic disinfectants and cleaning agents in the long term.

Never use aggressive or abrasive cleaning agents.



Do not use acetone and acetone-based cleaning agents to clean the device, as they could damage the surfaces.

-  Contaminated parts with which the patient has come into contact during the examination (forehead rest and face shield) should be cleaned with a disinfectant approved for the purpose. These parts are designed to be wiped down using mild cleaning agents and disinfectants such as suds, disinfectants based on quaternary ammonium compounds (0.2 %), glutaral (2 %) or isopropanol (60 %).

### **Cleaning the exterior components**

- All parts of the casing may be wiped off with a moist but not drip-wet cloth. Wipe off any marks or stains with distilled water, to which a drop of household washing up liquid has been added.

 Do not use any volatile cleansers, solvents or benzene, etc.

 Ensure that no moisture penetrates the device during cleaning and disinfection.

 Please use a soft, dry cloth to clean glass surfaces after having blown off any dust and other dirt.