

Practice social distancing with your patient exams.

Using ZEISS CIRRUS OCT



In this era of COVID-19, healthcare providers are required to take extra precautions to prevent the spread of pathogens by practicing social distancing.

ZEISS wants to support you with technologies and solutions that help you provide exceptional care to your patients, while also keeping everyone's health and safety in mind.


To facilitate a safe social distance during patient exams, ZEISS CIRRUS users have the option to use peripheral devices and accessories to extend the distance between operator and patient using the suggested exam configurations.



Seeing beyond

Connecting your device to an external monitor

If you would like to operate your CIRRUS at a safe social distance, you have the option to use peripheral accessories to help you do so. All you need is an external monitor, keyboard, and mouse.



This option is intended for the Operator and Patient to be in the same room at a safe distance.

Recommended Cable Lengths

- Passive VGA cables (five meters/15 feet or less is recommended)
- Passive DVI cables (dual-link cables, 15 meters/49 feet or less is recommended)
- Passive DisplayPort cables (five meters/15 feet or less is recommended)
- Passive HDMI cables (20 meters/65 feet or less is recommended)

Recommended Keyboard Extension

You may choose to extend the length of your corded keyboard cable using a passive USB extension cable. It's recommended that you limit the length of your extension cable to a maximum of six meters/19 feet to ensure reliability.

Wireless Keyboard and Mouse

You may also opt to forgo cables and use wireless accessories.

Ports and Compatibility

	CIRRUS 5000/500	CIRRUS 6000
Ports	VGA Port Display Port DVI Port	HDMI Port Display Port DVI Port
Video Signal	1280 x 1024; 60Hz; 4:3 aspect ratio	1920 x 1080; 60Hz (1080P); 16:9 aspect ratio
Standard Video Interface	DVI	HDMI or DVI
Dual Monitor Mode	Mirror	Mirror

For more information and additional resources visit our MED Support Now Website: www.zeiss.com/med-support-now