|                     |                                            |    |
|---------------------|--------------------------------------------|----|
| fr                  | NOTICE D'UTILISATION – PAINLESS PLUG®      | 2  |
| en                  | INSTRUCTION FOR USE – PAINLESS PLUG®       | 4  |
| en-US               | INSTRUCTION FOR USE – SNUG PLUG            | 6  |
| de                  | GEBRAUCHSANWEISUNG – PAINLESS PLUG®        | 8  |
| it                  | ISTRUZIONI PER L'USO – PAINLESS PLUG®      | 10 |
| es                  | INSTRUCCIONES DE USO – PAINLESS PLUG®      | 12 |
| pt                  | INSTRUÇÕES DE USO – PAINLESS PLUG®         | 14 |
| nl                  | GEBRUIKSAANWIJZING – PAINLESS PLUG®        | 16 |
| da                  | BRUGSVEJLEDNING – PAINLESS PLUG®           | 18 |
| fi                  | KÄYTTÖOHJEET – PAINLESS PLUG®              | 20 |
| sv                  | BRUKSANVISNING – PAINLESS PLUG®            | 22 |
| el                  | ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ – PAINLESS PLUG®            | 24 |
| et                  | KASUTUSJUHEND – PAINLESS PLUG®             | 26 |
| lv                  | LIETOSANAS INSTRUKCIJA – PAINLESS PLUG®    | 28 |
| lt                  | NAUDOJIMO INSTRUKCIJOS – PAINLESS PLUG®    | 30 |
| pl                  | INSTRUKCJA UŻYTKOWANIA – PAINLESS PLUG®    | 32 |
| sl                  | NAVODILA ZA UPORABO – PAINLESS PLUG®       | 34 |
| sk                  | NAVOD NA POUZIVANIE – PAINLESS PLUG®       | 36 |
| cs                  | NAVOD K POUZITI – PAINLESS PLUG®           | 38 |
| ro                  | INSTRUCȚIUNI DE UTILIZARE – PAINLESS PLUG® | 40 |
| bg                  | ИНСТРУКЦИИ ЗА УПОТРЕБА – PAINLESS PLUG®    | 42 |
| hr                  | UPUTE ZA UPOTREBU – PAINLESS PLUG®         | 44 |
| hu                  | HASZNALATI UTASITAS – PAINLESS PLUG®       | 46 |
| no                  | BRUKSANVISNING – PAINLESS PLUG®            | 48 |
| tr                  | KULLANMA TALIMATI – PAINLESS PLUG®         | 50 |
| zh-TW               | 使用說明 – PAINLESS PLUG®                      | 52 |
| pt-BR               | INSTRUÇÕES DE USO – PAINLESS PLUG®         | 54 |
| zh-CN               | 使用说明 - 泪点塞 (painless 式)                    | 56 |
| ar                  | تعليمات الاستخدام – PAINLESS PLUG®         | 58 |
| Illustrations ..... |                                            | 60 |
| Symbols .....       |                                            | 61 |
| IFU-2138-307-83_I   |                                            |    |

FCI S.A.S.  
20-22, rue Louis Armand  
75015 PARIS – FRANCE

☎ : + 33 (0) 1 53 98 98 98

📠 : + 33 (0) 1 53 98 98 99

@ : fci@fci.fr  
🌐 : www.fci.fr

Distributed in U.S.A. by:

FCI Ophthalmics Inc.  
30 Corporate Park Drive  
Suite # 310/320  
Pembroke, MA 02359

☎ : 800-932-4202

📠 : 781-826-9062

@ : info@fci-ophthalmics.com  
🌐 : www.fci-ophthalmics.com

R<sub>x</sub> Only

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## PAINLESS PLUG®

### Description and packaging

Product scope of this IFU: S2.4001 / S2.4002 : PAINLESS PLUG®

PAINLESS PLUG® is a silicone plug preloaded on a « Plug inserter » made of polypropylene and stainless steel.

The plug (Fig. A) is composed of:

- a bulb (1) to prevent the plug from spontaneous expulsion,
- a body (2) that connects the collarette to the bulb and includes a closed canal for plug placement,
- a collarette (3) that caps the punctal opening in order to prevent intra-canalicular migration of the plug.

The plug is provided preloaded and stretched on a « Plug inserter » (Fig. B), made of a metallic tip and a release mechanism that allows a progressive release by squeezing with the thumb and index finger.

“One size fits all” device. Dimensions of the released plug are presented in Fig. C.

The device is sterilized using ethylene oxide in its final simple blister packaging. A green indicator shows that the product followed a sterilization cycle validated by FCI.

The device is supplied by box of 1 or box of 2. The device is for single use only and is not resterilizable.

### Intended purpose

The device is designed to close the lacrimal puncta in the treatment of dry eye.

### Indications

Medical conditions: dry eye

Intended patient population: children (12 months and older) and adults

Body parts in contact with the device: lacrimal puncta, canaliculi, palpebral conjunctiva, eyelids

Intended users: professional users and trained ophthalmic surgeons

Intended use environment: consulting room

### Contraindications

Utilization of the device is contraindicated in the following cases:

- impermeability of the lacrimal ducts,
- abnormalities of the eyelids with absence of nocturnal palpebral occlusion,
- senile conjunctival hyperlaxity covering the inferior punctum,
- a positive Bengal Rose test beyond the area of the palpebral fissure,
- rupture of the punctum.

### Warnings

The device is for single use and must not be resterilized. Reuse and/or resterilization may compromise device performance, which could cause serious harm to the patient's health and safety.

Check packaging integrity before use, and do not use the device in the event of packaging being damaged, unintentionally opened before use; or if the packaging is exposed to environmental conditions outside of those specified.

Do not use the product if the indicator is not green, as it could endanger the patient. Do not use after the expiration date shown on the package.

The implantable part of the device is made entirely of medical grade silicone.

The device is usually left in place for 2 or 3 months.

However, the implantation duration of the device should be controlled and adapted to the healing of the patient, as observed by the physician.

The guaranteed maximum residual rate of ethylene oxide the device can release to the patient is 4 µg. The residual quantity of ethylene oxide corresponds to the maximum quantity of residue that may leach over the lifetime of use.

In case the user identifies faulty labelling information, it shall be relayed to FCI S.A.S. and the device should not be used if the patient's safety is threatened.

### Precautions for use

The device must be stored at temperatures between 5°C/41°F and 35°C/95°F and humidity between 30% and 75%, away from sunlight.

The device may be exposed to temperature variations between -15°C/5°F and 50°C/122°F and humidity variations between 15% and 90% during transport for 10 days.

**Adverse effects**

As in any type of surgery, there are risks linked to the material or to developments of the initial pathology.  
Potential complications associated with the implantation of the device include, but are not limited to the following:

- canaliculitis
- plug extrusion
- granuloma formation
- epiphora
- rupture of the punctum
- temporary pruritus
- inflammation
- canalicular stenosis
- conjunctival erosion
- infection
- plug loss
- conjunctival hyperemia
- irritation
- foreign body sensation
- discomfort

Any complication that occurs in relation with the device should be notified to FCI S.A.S. and the competent authority of the Member State in which the user and/or the patient is established, or, for users and/or patients not subject to EU requirements, to FCI S.A.S. and to the applicable regulatory authority in conformity with the requirements of that regulatory authority.

**Instructions for use**Insertion:

Prior dilation of the lacrimal punctum is recommended.

Insertion of the device using the « Plug inserter » (Fig.D):

Grasp eyelid with one hand to stabilize it horizontally and expose the punctal opening.

Insert the plug in the punctum until the plug inserter is abutted against the eyelid. (1)

Squeeze the inserter handles together until full contact with the body of the plug inserter is achieved. The gripper retracts and releases the plug which returns to its final shape. (2)

Keep squeezing the inserter before and while removing the plug inserter from the punctum. (3)

Removal:

Use blunt forceps to grasp plug below the collarette and gently pull it out from the lacrimal puncta.

**Performances, mode of action and clinical benefits**

The plug obstructs the lacrimal punctum to prevent the tears from draining away from the eye to the lacrimal canaliculi, in order to treat dry eye.

PAINLESS PLUG® is delivered in an extended position thereby elongating the plug to facilitate insertion.

The Summary of Safety and Clinical Performances is available in the European database on medical devices (Eudamed: <https://ec.europa.eu/tools/eudamed>), where it is linked to the Basic UDI-DI.

Basic UDI-DI: 3700773707079M

**Compatibility with other devices**

The implanted part of the device is MR safe.

**Information to be communicated to the patient**

There may be local discomfort for two hours after the procedure. Do not rub the eye to let the device take its place.

Do not rub the eye and do not touch the inner corner of the eye during the time of the treatment.

An implant card is provided in the box of the product. It is intended to be filled out by the healthcare providers and given to the patient. It allows the traceability of the implanted product. Patient information is made available on the internet. The website address is printed on the implant card.

Any serious incident that occurs in relation with the device should be notified to FCI S.A.S. and the competent authority of the State in which the patient is established.

The user must inform the patient of the contraindications and potential complications related to the device that are set out in the present document.

**Waste disposal**

| With household waste or recycle                         | With hospital waste (contaminated products) |
|---------------------------------------------------------|---------------------------------------------|
| Outer box, instructions for use and traceability labels | Plug                                        |
| Blister                                                 | Inserter                                    |

Use caution with parts of the devices which present the following risks: pricks, infection, microbial risk.

**En-US**

Instructions for use: 2138-307-83\_I  
Date of revision of the instructions for use : 2025-03

**SNUG PLUG****Description and packaging**

Product scope of this IFU: S2.4002U SNUG PLUG

SNUG PLUG is a silicone plug preloaded on a « Plug inserter » made of polypropylene and stainless steel.

The plug (Fig. A) is composed of:

- a bulb (1) to prevent the plug from spontaneous expulsion,
- a body (2) that connects the collarette to the bulb and includes a closed canal for plug placement,
- a collarette (3) that caps the punctal opening in order to prevent intra-canalicular migration of the plug.

The plug is provided preloaded and stretched on a « Plug inserter » (Fig. B), made of a metallic tip and a release mechanism that allows a progressive release by squeezing with the thumb and index finger.

“One size fits all” device. Dimensions of the released plug are presented in Fig. C.

The device is sterilized using ethylene oxide in its final simple blister packaging. A green indicator shows that the product followed a sterilization cycle validated by FCI.

The device is supplied by box of 1 or box of 2. The device is for single use only and is not resterilizable.

**Intended purpose**

The device is designed to close the lacrimal puncta in the treatment of dry eye.

**Indications**

Medical conditions: dry eye

Intended patient population: children (12 months and older) and adults

Body parts in contact with the device: lacrimal puncta, canaliculi, palpebral conjunctiva, eyelids

Intended users: professional users and trained ophthalmic surgeons

Intended use environment: consulting room

**Contraindications**

Utilization of the device is contraindicated in the following cases:

- impermeability of the lacrimal ducts,
- abnormalities of the eyelids with absence of nocturnal palpebral occlusion,
- senile conjunctival hyperlaxity covering the inferior punctum,
- a positive Bengal Rose test beyond the area of the palpebral fissure,
- rupture of the punctum

**Warnings**

The device is for single use and must not be resterilized. Reuse and/or resterilization may compromise device performance, which could cause serious harm to the patient's health and safety.

Check packaging integrity before use, and do not use the device in the event of packaging being damaged, unintentionally opened before use; or if the packaging is exposed to environmental conditions outside of those specified.

Do not use the product if the indicator is not green, as it could endanger the patient. Do not use after the expiration date shown on the package.

The implantable part of the device is made entirely of medical grade silicone.

The device is usually left in place for 2 or 3 months.

However, the implantation duration of the device should be controlled and adapted to the healing of the patient, as observed by the physician.

The guaranteed maximum residual rate of ethylene oxide the device can release to the patient is 4 µg. The residual quantity of ethylene oxide corresponds to the maximum quantity of residue that may leach over the lifetime of use.

In case the user identifies faulty labelling information, it shall be relayed to FCI S.A.S. and the device should not be used if the patient's safety is threatened.

**Precautions for use**

The device must be stored at temperatures between 5°C/41°F and 35°C/95°F and humidity between 30% and 75%, away from sunlight.

The device may be exposed to temperature variations between -15°C/5°F and 50°C/122°F and humidity variations between 15% and 90% during transport for 10 days.

**Adverse effects**

As in any type of surgery, there are risks linked to the material or to developments of the initial pathology.  
Potential complications associated with the implantation of the device include, but are not limited to the following:

- canaliculitis
- plug extrusion
- granuloma formation
- epiphora
- rupture of the punctum
- temporary pruritus
- inflammation
- canalicular stenosis
- conjunctival erosion
- infection
- plug loss
- conjunctival hyperemia
- irritation
- foreign body sensation
- discomfort

Any complication that occurs in relation with the device should be notified to FCI S.A.S. and the competent authority of the Member State in which the user and/or the patient is established, or, for users and/or patients not subject to EU requirements, to FCI S.A.S. and to the applicable regulatory authority in conformity with the requirements of that regulatory authority.

**Instructions for use**Insertion:

Prior dilation of the lacrimal punctum is recommended.

Insertion of the device using the « Plug inserter » (Fig.D):

Grasp eyelid with one hand to stabilize it horizontally and expose the punctal opening.

Insert the plug in the punctum until the plug inserter is abutted against the eyelid. (1)

Squeeze the inserter handles together until full contact with the body of the plug inserter is achieved. The gripper retracts and releases the plug which returns to its final shape. (2)

Keep squeezing the inserter before and while removing the plug inserter from the punctum. (3)

Removal:

Use blunt forceps to grasp plug below the collarette and gently pull it out from the lacrimal puncta.

**Performances, mode of action and clinical benefits**

The plug obstructs the lacrimal punctum to prevent the tears from draining away from the eye to the lacrimal canaliculi, in order to treat dry eye.

SNUG PLUG is delivered in an extended position thereby elongating the plug to facilitate insertion.

The Summary of Safety and Clinical Performances is available in the European database on medical devices (Eudamed: <https://ec.europa.eu/tools/eudamed>), where it is linked to the Basic UDI-DI.

Basic UDI-DI: 3700773707079M

**Compatibility with other devices**

The implanted part of the device is MR safe.

**Information to be communicated to the patient**

There may be local discomfort for two hours after the procedure. Do not rub the eye to let the device take its place.

Do not rub the eye and do not touch the inner corner of the eye during the time of the treatment.

An implant card is provided in the box of the product. It is intended to be filled out by the healthcare providers and given to the patient. It allows the traceability of the implanted product. Patient information is made available on the internet. The website address is printed on the implant card.

Any serious incident that occurs in relation with the device should be notified to FCI S.A.S. and the competent authority of the State in which the patient is established.

The user must inform the patient of the contraindications and potential complications related to the device that are set out in the present document.

**Waste disposal**

| With household waste or recycle                         | With hospital waste (contaminated products) |
|---------------------------------------------------------|---------------------------------------------|
| Outer box, instructions for use and traceability labels | Plug                                        |
| Blister                                                 | Inserter                                    |

Use caution with parts of the devices which present the following risks: pricks, infection, microbial risk.

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