

# Probability of surviving fast progression and eye growth reversal after 1-year of spectacle wear with cylindrical annular refractive elements

Poster #144 – A0444

Padmaja Sankaridurg<sup>1,2</sup>, Xiaoqin Chen<sup>4</sup>, Katharina Rifai<sup>1</sup>, Min Wu<sup>6</sup>, Cui Yu<sup>5</sup>, Arne Ohlendorf<sup>1</sup>, Christina Boeck-Maier<sup>1</sup>, Siegfried Wahl<sup>1,3</sup>, Youhua Yang<sup>7</sup>, Yi Zhu<sup>7</sup>, Lihua Li<sup>4</sup>

<sup>1</sup> ZEISS Vision Care, Carl Zeiss Vision International GmbH, Aalen, Germany; <sup>2</sup> School of Optometry and Vision Science, University of New South Wales, Sydney, NSW, Australia; <sup>3</sup> Institute for Ophthalmic Research, Eberhard Karls University of Tübingen, Tübingen, Germany; <sup>4</sup> Tianjin Eye Hospital Optometric Center, Tianjin, Tianjin, China; <sup>5</sup> Shenyang He Eye Hospital, Shenyang, Liaoning, China; <sup>6</sup> Beijing Tongren Vision Care, Beijing Tongren Hospital CMU, Beijing, Beijing, China; <sup>7</sup> ZEISS Vision Care, Carl Zeiss Vision (Guangzhou) Ltd., Guangzhou, Guangdong, China

## Purpose

Faster progression of myopia (annual change of  $\leq -0.50D$  or higher) is associated with greater risk of: visual impairment, more frequent optical device changes, more severe myopia and myopia-related complications in adult life.<sup>1-3</sup> We analyzed the likelihood of progression of  $-0.75D$  or worse after 1-year of lens wear and incidents of axial length (AL) reversal with single vision spectacles compared to MyoCare spectacle lenses incorporating cylindrical annular refractive elements (CARE).

## Methods

- Interim analysis of 12-month data from an ongoing 2-year prospective, double-masked, multi-center clinical trial (NCT05288335) (Figure 1)




Figure 1. Overview of the randomized controlled clinical trial.

- Change in spherical equivalent refractive error (SE) and axial length (AL) from baseline was determined.
- Survival analysis was determined for fast annual progression  $\leq -0.75D$  (i.e.,  $-0.75D$  or worse) using Kaplan-Meier Survival analysis with a log-rank test for differences between groups.
- Percent eyes exhibiting AL reversal at 6 and 12 months were determined.

## Results

- At baseline, there were no differences between the groups for any of the characteristics. Participant flow through the trial is presented in Figure 2.




Figure 2. Participant flow over 12 months.

Table 1. Baseline characteristics

Characteristic	MyoCare (n = 80)	MyoCare S (n = 80)	Single Vision (n = 80)
Age (years)	9.9 ± 1.7	9.8 ± 1.7	9.8 ± 1.6
Gender (M/F)	45(36.6%)/35(29.9%)	39(31.7%)/41(35.0%)	39(31.7%)/41(35.0%)
Parental Myopia			
None	14(20.5%)	21(31.3%)	18(26.9%)
One	26(38.2%)	23(34.3%)	23(34.3%)
Both	28(41.2%)	23(34.3%)	26(38.8%)
Cycloplegic SE(D) RE	-2.23 ± 0.98	-2.30 ± 1.06	-2.31 ± 1.01
AL (mm) RE	24.44 ± 0.73	24.34 ± 0.74	24.43 ± 0.73

## Probability of surviving fast progression ( $\leq -0.75D$ )

- Of the SV wearers, probability of surviving fast progression was 0.896 and 0.493 at 6 and 12 months (Figure 3).
- In comparison, both MyoCare lenses had a higher probability of survival that was significant compared to SV but similar between MyoCare and MyoCare S lenses at both 6 and 12 months.




Figure 3. Probability of surviving progression of  $\leq -0.75D$ .

## Reversal of axial length

- A small number of eyes demonstrated reversal at both 6 and 12 months. Reversal was greater with MyoCare lenses compared to SV (Figure 4).
- MyoCare > MyoCare S > SV but not significant.




Figure 4. Percent eyes demonstrating axial length reversal.

## Discussion

- Estimated annual myopia progression in Asian children with myopia (mean age 9.3 yrs) was  $-0.82D$  (95% CI,  $-0.71$  to  $-0.93D$ ).<sup>4</sup> In children wearing SV contact lenses, 41% had fast progression ( $-0.75D$  or worse).<sup>3</sup>
- In this trial, after 12 months of SV wear, probability of surviving fast progression was 0.493 or, one of every two SV wearing eyes experienced fast progression. In comparison, nearly 83% and 86% of MyoCare and MyoCare S wearing eyes survived fast progression. The differences in surviving fast progression were not different between MyoCare and MyoCare S, but significantly different compared to SV. Furthermore, at both 6 and 12 months, more eyes wearing MyoCare and MyoCare S exhibited axial length reversal.
- Our results indicate slower eye growth with MyoCare and MyoCare S compared to SV. AL reversal or shortening was reported with orthokeratology, atropine and low-level red-light therapy.<sup>5-7</sup> In the current study, AL shortening was greater at 6 months with a slight acceleration in eye growth thereafter. It has been suggested a thickened choroid, shrunken vitreous chamber or possibly a thinner cornea may be responsible.<sup>7</sup>

## Conclusion

Probability of surviving fast progression was greater with MyoCare than SV. Whilst 51% of SV eyes experienced fast progression ( $-0.75D$  or worse) in myopia after 12 months, only 17% of MyoCare and 14% of MyoCare S wearing eyes experienced fast progression.

## References

- Myles W, Dunlop C, McFadden SA. The Effect of Long-Term Low-Dose Atropine on Refractive Progression in Myopic Australian School Children. *J Clin Med*. 2021 Apr 1;10(7):1444.
- Shih YF, Chen CH, Chou AC, et al. Effects of different concentrations of atropine on controlling myopia in myopic children. *J Ocul Pharmacol Ther*. 1999 Feb;15(1):85-90.
- Sankaridurg P, Bakaraju RC, Naduvilath T, et al. Myopia control with novel central and peripheral plus contact lenses and extended depth of focus contact lenses: 2 year results from a randomised clinical trial. *Ophthalmic Physiol Opt*. 2019 Jul;39(4):294-307.
- Donovan L, Sankaridurg P, Ho A, Naduvilath T, Smith EL 3rd, Holden BA. Myopia progression rates in urban children wearing single-vision spectacles. *Optom Vis Sci*. 2012 Jan;89(1):27-32.
- Hu Y, Ding X, Jiang J et al. Long-Term Axial Length Shortening in Myopic Orthokeratology: Incident Probability, Time Course, and Influencing Factors. *Invest Ophthalmol Vis Sci*. 2023;64(15):37.
- Cho H, Seo Y, Han SH, Han J. Factors Related to Axial Length Elongation in Myopic Children Who Received 0.05% Atropine Treatment. *J Ocul Pharmacol Ther*. 2022 Dec;38(10):703-708.
- Wang W, Jiang Y, Zhu Z, Zhang S, et al. Axial Shortening in Myopic Children after Repeated Low-Level Red-Light Therapy: Post Hoc Analysis of a Randomized Trial. *Ophthalmol Ther*. 2023 Apr;12(2):1223-1237.

