

## **PR Fact Sheet**

## ZEISS VISUCORE 500 at a Glance

P	
What is ZEISS	ZEISS VISUCORE 500 is a highly precise, automated, combined and
VISUCORE 500?	thus space-saving objective-subjective refraction unit.
	Thanks to smart control modes and user-friendly operation, the
	system makes refraction faster and easier and it can be operated by
	all employees – regardless of their level of training and experience.
How much space	The ZEISS VISUCORE 500 unit occupies an area of only eight square
does the system	feet, respectively 0.8 square meters.
require?	
What are the	ZEISS VISUCORE 500 consists of ZEISS VISUSCREEN 100
component parts of	Compact and ZEISS VISUCORE 500 Refractor. An ergonomic table is
ZEISS VISUCORE	also an integral part of the entire system.
500?	
	The ZEISS VISUSCREEN 100 Compact vision testing unit contains
	various digital optotype charts (to test visual acuity and binocular
	vision). A distance of five meters is simulated in a distance of a mere
	one meter.
	The ZEISS VISUCORE 500 Refractor combines a binocular wavefront
	autorefractor, which enables ZEISS i.Scription, a keratometer for
	central corneal values and a digital phoropter with ZEISS-developed
	workflows (as well as red-green separation for binocular vision tests).
	Thanks to the integrated dual cross-cylinder lenses, both helical
	positions can be seen simultaneously during cylinder determination.
	This makes comparability easier for consumers as they do not have to
	remember any visual impression.
	A digital vision chart for testing near vision is also integrated – with a
	setting functionality for the appropriate distance of 25 to 70 cm placed
	directly on the device.
	ZEISS VISUCORE 500 operation is seamlessly integrated into the
	ZEISS VISUCORE 500 operation is seattlessly integrated into the ZEISS VISUCONSULT 500 app.



How can subjective	The mirrors built into ZEISS VISUSCREEN 100 Compact and a
refraction be	special lens make it possible to simulate the optically necessary
correct in such a	distance of five meters.
small space?	
What control modes	Three modes are available – with 1 and 2 being familiar to most eye
are integrated into	care professionals from previous ZEISS refraction applications:
the software? Where is ZEISS VISUCORE 500	<ol> <li>Workflow Mode: This comes with a predefined set of workflows and individual refraction modes can be saved here. Settings on the instrument are made automatically – depending on your selection of a test you previously saved.</li> <li>Freestyle Mode: This professional mode is used in special situations. All available charts can be selected individually. Any element of the phoropter can be steered individually.</li> <li>The new Guided Mode: This enables objective and subjective refraction to be done even by relatively inexperienced personnel. The software guides the user in detail through the entire process; provides specific instructions on how to interact with the person undergoing testing; and indicates questions that need to be asked; as well as the corresponding answers and the action to be taken as a result of the answers.</li> <li>ZEISS VISUCORE 500 is available on the US market since spring 2023. Availability will be gradually extended to other markets in the</li> </ol>
available?	near future.

VISUSTORE, VISUSCREEN, VISUCONSULT and VISUCORE are either trademarks or registered trademarks of Carl Zeiss AG or any other companies of the ZEISS Group in Germany and other countries.

Status: July 2023