

EARN
\$15
ON
EVERY PAIR



GREAT LENS. GREAT REWARDS.

Managed Care has never been more rewarding.

With AO Easy, your managed care patients will get an outstanding new progressive design – and you'll earn

\$15 on every pair*.
(Including private pay jobs.)

Use those funds for equipment, marketing, motivating your staff – whatever you think will help your practice grow.



Earning your rewards is **EASY.**

- **Order AO Easy online***
Your SOLA eRewards account is credited automatically.
- **View your completed orders**
by logging on to www.SOLAEwards.com
- **Get your rewards**
Click *SUBMIT* and SOLA will send a check directly to you in 6 to 8 weeks.

For more information go to:
www.sola.com/aoeasy/

Do the math!

If you dispense...	With AO Easy, you can earn...
5 progressive pairs per week	\$ 3,900 more per year
10 progressive pairs per week	\$ 7,800 more per year
15 progressive pairs per week	\$11,700 more per year

Keep earning!

SOLA eRewards is an ongoing program.* Keep ordering AO Easy and you'll keep earning rewards.

*Includes both AO Easy insurance and private pay jobs ordered online at www.eyefinity.com. No Medicare/Medicaid. All orders must be submitted online via www.eyefinity.com to be eligible for SOLA eRewards. To receive payment, participants must redeem through www.SOLAEwards.com. SOLA reserves the right, at its discretion, to end the promotion without notice. Void where prohibited by law. For complete rules visit www.SOLAEwards.com. Offer subject to change without notice.



GREAT LENS. GREAT REWARDS.

ao easy.

easy life.

easy Wear


- Smooth distribution of power and astigmatism assures wearing comfort.
- Balanced and unified viewing zones.

easy Adapt

- The AO Easy design significantly reduces the peripheral swim and distortion that can make adaptation difficult.
- The smooth transitions assure a natural vision experience – even for first-time progressive wearers.

easy Choice

- Optimized corridor length fits a greater range of frames (18mm minimum fitting height).
- Available in a full range of premium materials, and with Transitions®.



Availability

Material	Base Curves	Rx Range	Add Powers
AO Easy HD 1.67 High Index	Available Rx Teflon® coated only	-12.00 to +8.00 D	1.00 to 3.50
AO Easy HD 1.67 Transitions® V Gray	Available Rx Teflon® coated only	-12.00 to +8.00 D	1.00 to 3.50
NEW 1.67 High Index	1.00, 3.00, 5.00, 6.00, 7.00	-10.00 to +6.00 D	0.75 to 3.50
NEW 1.67 Transitions® V Gray	1.00, 3.00, 5.00, 6.00, 7.00	-10.00 to +6.00 D	0.75 to 3.50
1.6 High Index	1.00, 3.00, 5.00, 6.00, 7.00	-10.00 to +6.00 D	0.75 to 3.50
Polycarbonate Transitions® V Gray	1.50, 3.50, 5.50, 7.50	-9.00 to +5.00 D	1.00 to 3.00
Polycarbonate	1.50, 3.50, 5.50, 7.50	-9.00 to +5.00 D	1.00 to 3.00
Hard Resin Transitions® Brown	1.50, 3.50, 5.50, 7.50	-9.00 to +5.00 D	0.75 to 3.50
Hard Resin Transitions® Gray	1.50, 3.50, 5.50, 7.50	-9.00 to +5.00 D	0.75 to 3.50
Hard Resin	1.50, 3.50, 5.50, 7.50	-9.00 to +5.00 D	0.75 to 3.50



Teflon®
CLEAR COAT LENSES

AO Easy lenses are available with Teflon® Clear Coat.

Available from



www.sola.com