

# ZEISS EQ Workplace

Streamline your refractive cataract workflow.



Seeing beyond



EQ Workplace® from ZEISS allows you to calculate and select IOLs based on a comprehensive review of all relevant diagnostic data – anytime, anywhere. It can also be connected directly to your CALLISTO eye® from ZEISS to prepare your surgical assistance functions.

## Benefit from a new level of protection against never-events by streamlining and automating your processes:

- Biometric data from the IOLMaster® from ZEISS are automatically transferred and relevant data fields populate to calculate and select the IOL you need.
- You can easily verify whether the IOL you have selected pre-surgically is being implanted by comparing the IOL information with the data on the ZEISS CALLISTO eye screen.
- With ZEISS EQ Mobile app for iPad and iPhone, your biometric and surgical planning data is available anytime, anywhere.
- Personalize your IOL Constants



Scan  
to learn  
more.