

Which printers can I connect to a Matrix 800?

ZEISS Quick Help: Humphrey® Matrix 800

The Matrix 800 is a closed Linux operating system to which printer drivers cannot be added. Therefore, Postscript drivers are internally built into the software.

The Matrix 800 can utilize a peer-to-peer connection to Postscript printers using a network cable and crossover adapter or a crossover cable, which can be ordered by contacting ZEISS Customer Care at 1-800-341-6968. The instrument can also be connected to network printers either shared or directly.

Refer to the Humphrey Matrix user manual *Instructions for Use* for safe and effective operation of the instrument.



Carl Zeiss Meditec Inc.
5160 Hacienda Drive
Dublin, CA 94568
USA
www.zeiss.com/med



Carl Zeiss Meditec AG
Goeschwitzer Str. 51-52
07745 Jena
Germany
www.zeiss.com/med



SER.9254 Printed in the United States. CZ-IX/2017 United States edition: Only for sales in selected countries. The contents of the reference guide may differ from the current status of approval of the product or service offering in your country. Please contact our regional representatives for more information. Subject to change in design and scope of delivery and as a result of ongoing technical development. Humphrey is either a trademark or registered trademark Carl Zeiss Meditec AG or other companies of the ZEISS Group in Germany and/or other countries. © Carl Zeiss Meditec, Inc., 2017. All rights reserved.