

## **synchrony EASY ADAPT (One) - Fact Sheets**

### **Table of Contents**

#### **EASY ADAPT (One)**

|                                                                       |                |
|-----------------------------------------------------------------------|----------------|
| <i><b>EASY ADAPT (One) 1.67 High Index</b></i>                        | <b>1 – 3</b>   |
| <i><b>EASY ADAPT (One) 1.67 High Index PhotoFusion Gray/Brown</b></i> | <b>4 – 6</b>   |
| <i><b>EASY ADAPT (One) 1.59 Polycarbonate</b></i>                     | <b>7 – 9</b>   |
| <i><b>EASY ADAPT (One) 1.59 PhotoFusion Gray/Brown</b></i>            | <b>10 – 12</b> |
| <i><b>EASY ADAPT (One) 1.5</b></i>                                    | <b>13 – 15</b> |
| <i><b>EASY ADAPT (One) 1.5 PhotoFusion Gray/Brown</b></i>             | <b>16 – 18</b> |

# PRODUCT FACT SHEET



**Product Name** synchrony Semi-Finished Progressive Easy Adapt (One) 1.67

## 1. Customer Product Codes

| CAT Code: | LGC:  | LGC Description:                  | Notes:                              | FDA:                            |
|-----------|-------|-----------------------------------|-------------------------------------|---------------------------------|
| N85       | L6735 | SYNC SF P EAS A (ONE) 167 HC -B-  | Country of Origin required on Label | Yes, US FDA disclaimer Required |
| N84       | L6734 | SYNC SF P EAS A (ONE) 167 UNC -B- |                                     |                                 |

## 2. Product Range:

| Nominal Diameter (mm) | Nominal Base Curve (D) | Add Power Range (D) | Rx Range Recommended (D) |       |
|-----------------------|------------------------|---------------------|--------------------------|-------|
| 73                    | 1.25                   | 1.00 to 3.50        | -9.00 to                 | -5.25 |
| 78                    | 2.75                   | 1.00 to 3.50        | -5.00 to                 | -2.00 |
| 78                    | 4.00                   | 1.00 to 3.50        | -1.75 to                 | +1.25 |
| 78                    | 5.75                   | 1.00 to 3.50        | +1.50 to                 | +3.75 |
| 73                    | 6.50                   | 1.00 to 3.50        | +4.00 to                 | +6.00 |
| 73                    | 7.50                   | 1.00 to 3.50        | +6.25 to                 | +8.00 |

## 3. Material & Coating Specifications:

|                     |                                                  |
|---------------------|--------------------------------------------------|
| Type                | 1.67 MR7                                         |
| Refractive Index    | $n_d = 1.659$ $n_e = 1.664$                      |
| Abbe Number         | $v_d = 32$ $v_e = 32$                            |
| Density             | 1.35 g/cm <sup>3</sup>                           |
| UV Protection (ISO) | UVA : 100 %, UVB : 100 % 2.0 mm Centre Thickness |
| Coating Type        | Uncoated & Hard Coated                           |

## 4. Lens Geometry:

| Nominal Base | True Front Curve (D) (1.530) | Radii (mm) | Nominal Back Curve (D) | Centre Thickness (mm) ( $\pm 1.0$ mm) | Nominal Edge Thickness at 'A' (mm) |
|--------------|------------------------------|------------|------------------------|---------------------------------------|------------------------------------|
| 1.25         | 1.24                         | 427.42     | -6.00                  | 8.5                                   | 15                                 |
| 2.75         | 2.74                         | 193.43     | -6.00                  | 7.5                                   | 12.5                               |
| 4.00         | 4.01                         | 132.17     | -6.00                  | 6.5                                   | 10                                 |
| 5.75         | 5.77                         | 91.85      | -6.00                  | 7.5                                   | 8                                  |
| 6.50         | 6.54                         | 81.04      | -6.00                  | 8.5                                   | 8                                  |
| 7.50         | 7.54                         | 70.29      | -6.00                  | 9.5                                   | 8                                  |

Document Number: 00581.FS1

For External Distribution

Reviewed by: Geoff Riley, 10 Oct 2016

Page 1 of 3

# PRODUCT FACT SHEET



## Product Name



synchrony Semi-Finished Progressive Easy Adapt (One) 1.67

### 5. Surfacing: Recommended centre & edge thicknesses for all regions:

| Power Range (D)     | USA Market Minimum Centre Thickness | Int'l Market Minimum Centre Thickness (mm) | All Markets Minimum Edge Thickness (mm) |
|---------------------|-------------------------------------|--------------------------------------------|-----------------------------------------|
| High Minus to -2.00 | 1.5                                 | 1.2                                        |                                         |
| -1.75 to +2.00      | 2.0 or *                            | 2.0 or *                                   | 1.3 *                                   |
| +2.25 to High Plus  |                                     |                                            | 1.3                                     |

\* Thicknesses are based on whichever is the limiting factor

### 6. Product Markings:

|                                                                                                                                                                                                                                                                                     |                                                                                                                                                           |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Ink Markings:</b> <div style="display: flex; justify-content: space-around; align-items: center;">   </div> |                                                                                                                                                           |
| <b>Front Surface: Coding</b> <div style="text-align: center;">  </div>                                                                                                                           | <b>Back Surface: Coding</b> <div style="text-align: center;">  </div> |

Document Number: 00581.FS1

For External Distribution

Reviewed by: Geoff Riley, 10 Oct 2016

Page 2 of 3


# PRODUCT FACT SHEET



**Product Name**

synchrony Semi-Finished Progressive Easy Adapt (One) 1.67

## 7. Product Range (Minus Cylinder):



Document Number: 00581.FS1

For External Distribution

Reviewed by: Geoff Riley, 10 Oct 2016

Page 3 of 3

# PRODUCT FACT SHEET



## Product Name

synchrony Semi-Finished Progressive Easy Adapt (One) 1.67 Photochromic

### 1. Customer Product Codes

| CAT Code: | LGC:  | LGC Description:                        | Notes:                              | FDA:                            |
|-----------|-------|-----------------------------------------|-------------------------------------|---------------------------------|
| N86       | L6736 | SYNC SF P EAS A (ONE) 167 PH GRY HC AB- | Country of Origin required on Label | Yes, US FDA disclaimer Required |
| Y66       | L6797 | SYNC SF P EAS A (ONE) 167 PH BRN HC AB- |                                     |                                 |

### 2. Product Range:

| Nominal Diameter (mm) | Nominal Base Curve (D) | Add Power Range (D) | Rx Range Recommended (D) |       |
|-----------------------|------------------------|---------------------|--------------------------|-------|
| 73                    | 1.25                   | 1.00 to 3.50        | -9.00 to                 | -5.25 |
| 78                    | 2.75                   | 1.00 to 3.50        | -5.00 to                 | -2.00 |
| 78                    | 4.00                   | 1.00 to 3.50        | -1.75 to                 | +1.25 |
| 78                    | 5.75                   | 1.00 to 3.50        | +1.50 to                 | +3.75 |
| 73                    | 6.50                   | 1.00 to 3.50        | +4.00 to                 | +6.00 |
| 73                    | 7.50                   | 1.00 to 3.50        | +6.25 to                 | +8.00 |

### 3. Material & Coating Specifications:

|                     |                                                  |
|---------------------|--------------------------------------------------|
| Type                | <b>1.67 MR7</b>                                  |
| Refractive Index    | $n_d = 1.659$ $n_e = 1.664$                      |
| Abbe Number         | $v_d = 32$ $v_e = 32$                            |
| Density             | 1.35 g/cm <sup>3</sup>                           |
| UV Protection (ISO) | UVA : 100 %, UVB : 100 % 2.0 mm Centre Thickness |
| Coating Type        | <b>Hard Coated</b>                               |

### 4. Lens Geometry:

| Nominal Base | True Front Curve (D) (1.530) | Radii (mm) | Nominal Back Curve (D) | Centre Thickness (mm) ( $\pm 1.0\text{mm}$ ) | Nominal Edge Thickness at 'A' (mm) |
|--------------|------------------------------|------------|------------------------|----------------------------------------------|------------------------------------|
| 1.25         | 1.16                         | 456.89     | -6.00                  | <b>8.5</b>                                   | 15                                 |
| 2.75         | 2.66                         | 199.25     | -6.00                  | <b>7.5</b>                                   | 12.5                               |
| 4.00         | 3.92                         | 135.20     | -6.00                  | <b>6.5</b>                                   | 10                                 |
| 5.75         | 5.68                         | 93.31      | -6.00                  | <b>7.5</b>                                   | 8                                  |
| 6.50         | 6.45                         | 82.17      | -6.00                  | <b>8.5</b>                                   | 8                                  |
| 7.50         | 7.46                         | 71.04      | -6.00                  | <b>9.5</b>                                   | 8                                  |

Document Number: 00582.FS1

For External Distribution

Reviewed by: Geoff Riley, 22 Nov 2016

Page 1 of 3

# PRODUCT FACT SHEET



## Product Name




synchrony Semi-Finished Progressive Easy Adapt (One) 1.67 Photochromic

### 5. Surfacing: Recommended centre & edge thicknesses for all regions:

| Power Range (D)     | USA Market Minimum Centre Thickness | Int'l Market Minimum Centre Thickness (mm) | All Markets Minimum Edge Thickness (mm) |
|---------------------|-------------------------------------|--------------------------------------------|-----------------------------------------|
| High Minus to -2.00 | 1.5                                 | 1.2                                        |                                         |
| -1.75 to +2.00      | 2.0 or *                            | 2.0 or *                                   | 1.3 *                                   |
| +2.25 to High Plus  |                                     |                                            | 1.3                                     |

\* Thicknesses are based on whichever is the limiting factor

### 6. Product Markings:

|                                                                                                                                                                                                                                                                                     |                                                                                                                                                           |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Ink Markings:</b> <div style="display: flex; justify-content: space-around; align-items: center;">   </div> |                                                                                                                                                           |
| <b>Front Surface: Coding</b> <div style="text-align: center;">  </div>                                                                                                                           | <b>Back Surface: Coding</b> <div style="text-align: center;">  </div> |

Document Number: 00582.FS1

For External Distribution

Reviewed by: Geoff Riley, 22 Nov 2016

Page 2 of 3


# PRODUCT FACT SHEET



**Product Name**

synchrony Semi-Finished Progressive Easy Adapt (One) 1.67 Photochromic

## 7. Product Range (Minus Cylinder):



Document Number: 00582.FS1

For External Distribution

Reviewed by: Geoff Riley, 22 Nov 2016

Page 3 of 3

# PRODUCT FACT SHEET



|                     |                                                           |
|---------------------|-----------------------------------------------------------|
| <b>Product Name</b> | synchrony Semi-Finished Progressive Easy Adapt (One) 1.59 |
|---------------------|-----------------------------------------------------------|

## 1. Customer Product Codes

| CAT Code: | LGC:  | LGC Description:                 | Notes:                              | FDA:                       |
|-----------|-------|----------------------------------|-------------------------------------|----------------------------|
| 270       | L2073 | SYNC SF P EAS A (ONE) 159 HC --- | Country of Origin required on Label | US FDA disclaimer Required |
| G36       | L5031 | SYNC SF P EAS A (ONE) 159 PG --- |                                     |                            |

## 2. Product Range:

| Nominal Diameter (mm) | Nominal Base Curve (D) | Add Power Range (D) | Rx Range Recommended (D) |
|-----------------------|------------------------|---------------------|--------------------------|
| 75                    | 1.25                   | 1.00 to 3.00        | -9.00 to -5.25           |
| 75                    | 3.25                   | 1.00 to 3.00        | -5.00 to -2.00           |
| 75                    | 4.50                   | 1.00 to 3.00        | -1.75 to +1.25           |
| 75                    | 6.25                   | 1.00 to 3.00        | +1.50 to +3.75           |
| 75                    | 7.00                   | 1.00 to 3.00        | +4.00 to +6.00           |

## 3. Material & Coating Specifications:

|                            |                               |
|----------------------------|-------------------------------|
| <b>Type</b>                | <b>1.59 Polycarbonate</b>     |
| <b>Refractive Index</b>    | $n_d = 1.586$ $n_e = 1.590$   |
| <b>Abbe Number</b>         | $v_d = 30$ $v_e = 29$         |
| <b>Density</b>             | 1.20 g/cm <sup>3</sup>        |
| <b>UV Protection (ISO)</b> | UVA : 100 %,      UVB : 100 % |
| <b>Coating Type</b>        | <b>Hard Coated</b>            |

## 4. Lens Geometry:

| Nominal Base | True Front Curve (D) (1.530) | Radii (mm) | Nominal Back Curve (D) | Centre Thickness (± 1.0mm) | Nominal Edge Thickness at 'A' (mm) | Nominal Packaged Weight (g) |
|--------------|------------------------------|------------|------------------------|----------------------------|------------------------------------|-----------------------------|
| 1.25         | 1.28                         | 414.06     | -5.75                  | 11.6                       | 17.8                               | 83                          |
| 3.25         | 3.22                         | 164.60     | -5.75                  | 10.6                       | 14.1                               | 70                          |
| 4.50         | 4.51                         | 117.52     | -5.75                  | 8.2                        | 10.0                               | 53                          |
| 6.25         | 6.22                         | 85.21      | -5.75                  | 10.4                       | 9.8                                | 58                          |
| 7.00         | 7.00                         | 75.71      | -8.00                  | 8.6                        | 10.2                               | 55                          |

|                              |                                       |
|------------------------------|---------------------------------------|
| Document Number: 0102274.FS6 | Reviewed by: Geoff Riley, 14 Sep 2016 |
| For External Distribution    | Page 1 of 3                           |

# PRODUCT FACT SHEET



## Product Name

synchrony Semi-Finished Progressive Easy Adapt (One) 1.59

### 5. Surfacing: Recommended centre & edge thicknesses for all regions:

| Power Range (D)     | USA Market Minimum Centre Thickness | Int'l Market Minimum Centre Thickness (mm) | All Markets Minimum Edge Thickness (mm) |
|---------------------|-------------------------------------|--------------------------------------------|-----------------------------------------|
| High Minus to -2.00 | 2.0                                 | 2.0                                        |                                         |
| -1.75 to +2.00      | 2.0 or *                            | 2.0 or *                                   | 1.5 *                                   |
| +2.25 to High Plus  |                                     |                                            | 1.5                                     |

\* Thicknesses are based on whichever is the limiting factor

### 6. Product Markings:

|                              |                             |
|------------------------------|-----------------------------|
| <b>Ink Markings:</b>         |                             |
| <b>Front Surface: Coding</b> | <b>Back Surface: Coding</b> |

Document Number: 0102274.FS6

For External Distribution

Reviewed by: Geoff Riley, 14 Sep 2016

Page 2 of 3


# PRODUCT FACT SHEET



**Product Name**

synchrony Semi-Finished Progressive Easy Adapt (One) 1.59

## 7. Product Range (Minus Cylinder):



Document Number: 0102274.FS6

For External Distribution

Reviewed by: Geoff Riley, 14 Sep 2016

Page 3 of 3

# PRODUCT FACT SHEET



## Product Name

synchrony Semi-Finished Progressive Easy Adapt (One) 1.59  
Photochromic

## 1. Customer Product Codes

| CAT Code: | LGC:  | LGC Description:                     | Notes:                              | FDA:                       |
|-----------|-------|--------------------------------------|-------------------------------------|----------------------------|
| F37       | L6086 | SYNC SF P EAS A (ONE) 159 PH GRY --- | Country of Origin required on Label | US FDA disclaimer Required |
| F38       | L6087 | SYNC SF P EAS A (ONE) 159 PH BRN --- |                                     |                            |

## 2. Product Range:

| Nominal Diameter (mm) | Nominal Base Curve (D) | Add Power Range (D) | Rx Range Recommended (D) |
|-----------------------|------------------------|---------------------|--------------------------|
| 72                    | 1.25                   | 1.00 to 3.00        | -9.00 to -5.25           |
| 72                    | 3.25                   | 1.00 to 3.00        | -5.00 to -2.00           |
| 72                    | 4.50                   | 1.00 to 3.00        | -1.75 to +1.25           |
| 72                    | 6.25                   | 1.00 to 3.00        | +1.50 to +3.75           |
| 72                    | 7.00                   | 1.00 to 3.00        | +4.00 to +6.00           |

## 3. Material & Coating Specifications:

| Type                | 1.59 Polycarbonate            |
|---------------------|-------------------------------|
| Refractive Index    | $n_d = 1.586$ $n_e = 1.590$   |
| Abbe Number         | $v_d = 30$ $v_e = 29$         |
| Density             | 1.20 g/cm <sup>3</sup>        |
| UV Protection (ISO) | UVA : 100 %,      UVB : 100 % |
| Coating Type        | Uncoated & Hard Coated        |

## 4. Lens Geometry:

| Nominal Base | True Front Curve (D) (1.530) | Radii (mm) | Nominal Back Curve (D) | Centre Thickness (± 1.0mm) | Nominal Edge Thickness at 'A' (mm) |
|--------------|------------------------------|------------|------------------------|----------------------------|------------------------------------|
| 1.25         | 1.28                         | 414.06     | -5.75                  | 11.6                       | 17.3                               |
| 3.25         | 3.26                         | 162.58     | -5.75                  | 10.6                       | 13.8                               |
| 4.50         | 4.49                         | 118.04     | -5.75                  | 8.2                        | 9.8                                |
| 6.25         | 6.25                         | 84.80      | -5.75                  | 10.4                       | 9.8                                |
| 7.00         | 7.01                         | 75.61      | -8.00                  | 8.9                        | 10.8                               |

Document Number: 00236.FS4

For External Distribution

Reviewed by: Geoff Riley, 20 Jul 2016

Page 1 of 3

# PRODUCT FACT SHEET



## Product Name

synchrony Semi-Finished Progressive Easy Adapt (One) 1.59  
Photochromic

### 5. Surfacing: Recommended centre & edge thicknesses for all regions:

| Power Range (D)     | USA Market Minimum Centre Thickness | Int'l Market Minimum Centre Thickness (mm) | All Markets Minimum Edge Thickness (mm) |
|---------------------|-------------------------------------|--------------------------------------------|-----------------------------------------|
| High Minus to -2.00 | 2.0                                 | 2.0                                        |                                         |
| -1.75 to +2.00      | 2.0 or *                            | 2.0 or *                                   | 1.5 *                                   |
| +2.25 to High Plus  |                                     |                                            | 1.5                                     |

\* Thicknesses are based on whichever is the limiting factor

### 6. Product Markings:

| Ink Markings:         |                      |
|-----------------------|----------------------|
|                       |                      |
| Front Surface: Coding | Back Surface: Coding |
|                       |                      |

Document Number: 00236.FS4

For External Distribution

Reviewed by: Geoff Riley, 20 Jul 2016

Page 2 of 3


# PRODUCT FACT SHEET



**Product Name**

synchrony Semi-Finished Progressive Easy Adapt (One) 1.59  
Photochromic

## 7. Product Range (Minus Cylinder):



Document Number: 00236.FS4

For External Distribution

Reviewed by: Geoff Riley, 20 Jul 2016

Page 3 of 3

# PRODUCT FACT SHEET



|                     |                                                          |
|---------------------|----------------------------------------------------------|
| <b>Product Name</b> | synchrony Semi-Finished Progressive Easy Adapt (One) 1.5 |
|---------------------|----------------------------------------------------------|

## 1. Customer Product Codes

| CAT Code: | LGC:  | LGC Description:                 | Notes:                              | FDA:                       |
|-----------|-------|----------------------------------|-------------------------------------|----------------------------|
| 275       | L2075 | SYNC SF P EAS A (ONE) 15 HC ---  | Country of Origin required on Label | US FDA disclaimer Required |
| 375       | L2401 | SYNC SF P EAS A (ONE) 15 UNC --- |                                     |                            |

## 2. Product Range:

| Nominal Diameter (mm) | Nominal Base Curve (D) | Add Power Range (D) | Rx Range Recommended (D) |
|-----------------------|------------------------|---------------------|--------------------------|
| 80                    | 1.25                   | 1.00 to 3.50        | -9.00 to -5.25           |
| 80                    | 3.25                   | 1.00 to 3.50        | -5.00 to -2.00           |
| 80                    | 4.50                   | 1.00 to 3.50        | -1.75 to +1.25           |
| 80                    | 6.25                   | 1.00 to 3.50        | +1.50 to +3.75           |
| 80                    | 7.00                   | 1.00 to 3.50        | +4.00 to +6.00           |

## 3. Material & Coating Specifications:

|                            |                                   |
|----------------------------|-----------------------------------|
| <b>Type</b>                | <b>1.5 Polymer</b>                |
| <b>Refractive Index</b>    | $n_d = 1.499$ $n_e = 1.500$       |
| <b>Abbe Number</b>         | $v_d = 58$ $v_e = 56$             |
| <b>Density</b>             | 1.32 g/cm <sup>3</sup>            |
| <b>UV Protection (ISO)</b> | UVA : 94 %,      UVB : 100 %      |
| <b>Coating Type</b>        | <b>Uncoated &amp; Hard Coated</b> |

## 4. Lens Geometry:

| Nominal Base | True Front Curve (D) (1.530) | Radii (mm) | Nominal Back Curve (D) | Centre Thickness (mm) ( $\pm 1.0\text{mm}$ ) | Nominal Edge Thickness at 'A' (mm) |
|--------------|------------------------------|------------|------------------------|----------------------------------------------|------------------------------------|
| 1.25         | 1.24                         | 427.42     | -3.00                  | <b>13.0</b>                                  | 15.5                               |
| 3.25         | 3.29                         | 161.09     | -3.00                  | <b>15.0</b>                                  | 14.5                               |
| 4.50         | 4.54                         | 116.74     | -6.00                  | <b>9.0</b>                                   | 11.5                               |
| 6.25         | 6.32                         | 83.86      | -6.00                  | <b>9.5</b>                                   | 9.0                                |
| 7.00         | 7.08                         | 74.86      | -6.00                  | <b>11.0</b>                                  | 10.0                               |

Document Number: 5116515.FS1

For External Distribution

Reviewed by: Geoff Riley, 20 Jul 2016

Page 1 of 3

# PRODUCT FACT SHEET



|                     |                                                          |
|---------------------|----------------------------------------------------------|
| <b>Product Name</b> | synchrony Semi-Finished Progressive Easy Adapt (One) 1.5 |
|---------------------|----------------------------------------------------------|

## 5. Surfacing: Recommended centre & edge thicknesses for all regions:

| Power Range (D)     | USA Market Minimum Centre Thickness | Int'l Market Minimum Centre Thickness (mm) | All Markets Minimum Edge Thickness (mm) |
|---------------------|-------------------------------------|--------------------------------------------|-----------------------------------------|
| High Minus to -2.00 | 2.0                                 | 2.0                                        |                                         |
| -1.75 to +2.00      | 2.0 or *                            | 2.0 or *                                   | 1.5 *                                   |
| +2.25 to High Pus   |                                     |                                            | 1.5                                     |

\* Thicknesses are based on whichever is the limiting factor

## 6. Product Markings:

|                              |                             |
|------------------------------|-----------------------------|
| <b>Ink Markings:</b>         |                             |
|                              |                             |
| <b>Front Surface: Coding</b> | <b>Back Surface: Coding</b> |
|                              |                             |

|                              |                                       |
|------------------------------|---------------------------------------|
| Document Number: 5116515.FS1 | Reviewed by: Geoff Riley, 20 Jul 2016 |
| For External Distribution    | Page 2 of 3                           |


# PRODUCT FACT SHEET



**Product Name**

synchrony Semi-Finished Progressive Easy Adapt (One) 1.5

## 7. Product Range (Minus Cylinder):



Document Number: 5116515.FS1

For External Distribution

Reviewed by: Geoff Riley, 20 Jul 2016

Page 3 of 3

# PRODUCT FACT SHEET



## Product Name

synchrony Semi-Finished Progressive Easy Adapt (One) 1.5  
Photochromic

### 1. Customer Product Codes

| CAT Code: | LGC:  | LGC Description:                       | Notes:                              | FDA:                       |
|-----------|-------|----------------------------------------|-------------------------------------|----------------------------|
| Y63       | L6793 | SYNC SF P EAS A (ONE) 15 PH GRY HC A-- | Country of Origin required on Label | US FDA disclaimer Required |
| Y69       | L6800 | SYNC SF P EAS A (ONE) 15 PH BRN HC A-- |                                     |                            |

### 2. Product Range:

| Nominal Diameter (mm) | Nominal Base Curve (D) | Add Power Range (D) | Rx Range Recommended (D) |       |
|-----------------------|------------------------|---------------------|--------------------------|-------|
| 78                    | 1.25                   | 1.00 to 3.50        | -9.00 to                 | -5.25 |
| 78                    | 3.25                   | 1.00 to 3.50        | -5.00 to                 | -2.00 |
| 78                    | 4.50                   | 1.00 to 3.50        | -1.75 to                 | +1.25 |
| 78                    | 6.25                   | 1.00 to 3.50        | +1.50 to                 | +3.75 |
| 78                    | 7.00                   | 1.00 to 3.50        | +4.00 to                 | +6.00 |

### 3. Material & Coating Specifications:

| Type                | 1.5 Polymer                                                          |               |
|---------------------|----------------------------------------------------------------------|---------------|
| Refractive Index    | $n_d = 1.499$                                                        | $n_e = 1.500$ |
| Abbe Number         | $v_d = 58$                                                           | $v_e = 56$    |
| Density             | 1.32 g/cm <sup>3</sup>                                               |               |
| UV Protection (ISO) | Coated, Activated: 2.0mm Centre Thickness<br>UVA : 100%, UVB : 100 % |               |
| Coating Type        | Hard Coated                                                          |               |

### 4. Lens Geometry:

| Nominal Base | True Front Curve (D) (1.530) | Radii (mm) | Nominal Back Curve (D) | Centre Thickness (mm) ( $\pm 1.0$ mm) | Nominal Edge Thickness at 'A' (mm) |
|--------------|------------------------------|------------|------------------------|---------------------------------------|------------------------------------|
| 1.25         | 1.24                         | 427.42     | -3.00                  | 13.0                                  | 15.5                               |
| 3.25         | 3.27                         | 162.08     | -3.00                  | 15.0                                  | 14.5                               |
| 4.50         | 4.54                         | 116.74     | -6.00                  | 9.0                                   | 11.5                               |
| 6.25         | 6.35                         | 83.47      | -6.00                  | 9.5                                   | 9.0                                |
| 7.00         | 7.10                         | 74.65      | -6.00                  | 11.0                                  | 10.0                               |

Document Number: 00574.FS1

For External Distribution

Reviewed by: Geoff Riley, 13 Oct 2016

Page 1 of 3

# PRODUCT FACT SHEET



## Product Name



synchrony Semi-Finished Progressive Easy Adapt (One) 1.5  
Photochromic

### 5. Surfacing: Recommended centre & edge thicknesses for all regions:

| Power Range (D)     | USA Market Minimum Centre Thickness | Int'l Market Minimum Centre Thickness (mm) | All Markets Minimum Edge Thickness (mm) |
|---------------------|-------------------------------------|--------------------------------------------|-----------------------------------------|
| High Minus to -2.00 | 2.0                                 | 2.0                                        |                                         |
| -1.75 to +2.00      | 2.0 or *                            | 2.0 or *                                   | 1.5 *                                   |
| +2.25 to High Plus  |                                     |                                            | 1.5                                     |

\* Thicknesses are based on whichever is the limiting factor

### 6. Product Markings:

|                                                                                                                  |                                                                                                                  |
|------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------|
| <b>Ink Markings:</b>          |                                                                                                                  |
| <b>Front Surface: Coding</b>  | <b>Back Surface: Coding</b>  |

Document Number: 00574.FS1

For External Distribution

Reviewed by: Geoff Riley, 13 Oct 2016

Page 2 of 3


# PRODUCT FACT SHEET



## Product Name

synchrony Semi-Finished Progressive Easy Adapt (One) 1.5  
Photochromic

## 7. Product Range (Minus Cylinder):



Document Number: 00574.FS1

For External Distribution

Reviewed by: Geoff Riley, 13 Oct 2016

Page 3 of 3