








IOL-Injector Correlation Chart

IOL View















	IOL	Diopter Range & Injector / Cartridge	
4-Haptic MICS	CT SPHERIS® 209M		0.0 to +32.0 D AT.Shooter A1-2000 injector / VISCOJECT™-BIO 1.8 Cartridge Set or VISCOJECT-BIO 1.8 Injector Set or A6 injector Set
	CT ASPHINA® 409M		
	CT ASPHINA 509M		
	AT LISA® 809M		
	AT TORBI® 709M		-10.0 to +28.0 D sphere ; +1.0 to +4.0 D cylinder AT.Shooter A1-2000 injector / VISCOJECT™-BIO 1.8 Cartridge Set or VISCOJECT-BIO 1.8 Injector Set or A6 injector Set
			-10.0 to +28.0 D sphere ; +4.5 to +12.0 D cylinder VISCOJECT-BIO 2.2 Injector Set or AT.Shooter A1-2000 injector / VISCOJECT-BIO 2.2 Cartridge Set +28.5 to +32.0 D sphere ; +1.0 to +12.0 D cylinder VISCOJECT-BIO 2.2 Injector Set or AT.Shooter A1-2000 injector / VISCOJECT-BIO 2.2 Cartridge Set
	AT TORBI® 719M		-8.0 to +28.0 D spherical equivalent (SE)* ; +1.0 to +4.0 D cylinder AT.Shooter A1-2000 injector / VISCOJECT™-BIO 1.8 Cartridge Set or VISCOJECT-BIO 1.8 Injector Set or A6 injector Set
			-4.0 to +28.0 D spherical equivalent (SE)* ; +4.5 to +12.0 D cylinder VISCOJECT-BIO 2.2 Injector Set or AT.Shooter A1-2000 injector / VISCOJECT-BIO 2.2 Cartridge Set +28.5 to +32.0 D spherical equivalent (SE)* ; +1.0 to +12.0 D cylinder VISCOJECT-BIO 2.2 Injector Set or AT.Shooter A1-2000 injector / VISCOJECT-BIO 2.2 Cartridge Set
AT LISA toric 909M		-10.0 to +24.0 D sphere ; +1.0 to +12.0 D cylinder AT.Shooter A1-2000 injector / VISCOJECT™-BIO 1.8 Cartridge Set or VISCOJECT-BIO 1.8 Injector Set or A6 injector Set +24.5 to +32.0 D sphere ; +1.0 to +12.0 D cylinder VISCOJECT-BIO 2.2 Injector Set or AT.Shooter A1-2000 injector / VISCOJECT-BIO 2.2 Cartridge Set	

*Further preselected SE/cylinder combinations are available above and below the stated SE range.

IOL-Injector Correlation Chart

IOL View



	IOL	Diopter Range & Injector / Cartridge	
4-Haptic MICS	AT LISA tri toric 939M		-10.0 to +28.0 D sphere; +4.5 to +12.0 D cylinder and +28.5 to +32.0 D sphere ; +1.0 to +12.0 D cylinder VISCOJECT-BIO 2.2 Injector Set or AT.Shooter A1-2000 injector / VISCOJECT-BIO 2.2 Cartridge Set
	AT LISA tri toric 949M		-5.0 to +32.0 D spherical equivalent (SE); +4.5 to +12.0 D cylinder and +32.5 to +35.0 D spherical equivalent (SE) ; +1.0 to +12.0 D cylinder VISCOJECT-BIO 2.2 Injector Set or AT.Shooter A1-2000 injector / VISCOJECT-BIO 2.2 Cartridge Set
	AT LARA toric 929M		-4.0 to +32.0 D spherical equivalent (SE)*; +4.5 to +12.0 D cylinder VISCOJECT-BIO 2.2 Injector Set or AT.Shooter A1-2000 injector / VISCOJECT-BIO 2.2 Cartridge Set
	BLUEMIXS® Preloaded IOLs		
	CT ASPHINA 409MP		0.0 to +32.0 D BLUEMIXS 180 Injector
	CT ASPHINA 509MP		
	AT LISA 809MP		
	AT LISA tri 839MP		
	AT LARA 829MP		
	AT TORBI 709MP		-10.0 to +28.0 D sphere; +1.0 to +4.0 D cylinder BLUEMIXS 180 Injector
AT TORBI 719MP		-8.0 to +28.0 D spherical equivalent (SE)*; +1.0 to +4.0 D cylinder BLUEMIXS 180 Injector	
AT LISA tri toric 939MP		-10.0 to +28.0 D sphere; +1.0 to +4.0 D cylinder BLUEMIXS 180 Injector	
AT LISA tri toric 949 MP		-5.0 to +32.0 D spherical equivalent; +1.0 to +4.0 D cylinder BLUEMIXS 180 Injector	

*Further preselected SE/cylinder combinations are available above and below the stated SE range.

IOL-Injector Correlation Chart

IOL View



	IOL		Diopter Range & Injector / Cartridge	
3-Haptic	AT LARA toric 929MP		-8.0 D to +32.0 D spherical equivalent (SE)*; +1.0 to +4.0 D cylinder BLUEMIXS 180 Injector	
	CT SPHERIS® 204		-10.0 to +30.0 D AT.Shooter A1-2000 / VISCOJECT-BIO 2.2 Cartridge Set or VISCOJECT-BIO 2.2 Injector Set	
	CT ASPHINA® 404		+31.0 to +45.0 D (CT SPHERIS 204) +31.0 to +42.0 D (CT ASPHINA 404) ACCUJECT™ 3.0 Injector Set	
	CT SPHERIS 203		+8.0 to +32.0 D Forceps	
	<i>SkyJet® Preloaded IOLs</i>			
	CT SPHERIS 203P		+8.0 to +26.0 D SkyJet 2.8 Injector +27.0 to +30.0 D Forceps	
	CT ASPHINA 603P		+10.0 to +30.0 D SkyJet 2.8 Injector	

IOL-Injector Correlation Chart

IOL View





	IOL	Diopter Range & Injector / Cartridge	
1-Piece C-Loop	<i>Preloaded hydrophobic IOLs</i>		
	CT LUCIA 621P/PY		0.0 to +24.0 D BLUESERT 2.2 Injector Set
			+24.5 to +30.0 D BLUESERT 2.4 Injector Set
			+30.5 to +34.0 D BLUESERT 2.6 Injector Set
	CT LUCIA 221P		0.0 to +24.0 D BLUESERT 2.2 Injector Set
			+24.5 to +30.0 D BLUESERT 2.4 Injector Set
+30.5 to +34.0 D BLUESERT 2.6 Injector Set			
AT ELANA 841P		0.0 to +24.0 D BLUESERT 2.2 Injector Set	
		+24.5 to +30.0 D BLUESERT 2.4 Injector Set	
		+30.5 to +34.0 D BLUESERT 2.6 Injector Set	

IOL-Injector Correlation Chart

IOL View



	IOL	Diopter Range & Injector / Cartridge	
3 Piece	CT LUCIA 202		+4.0 to +34.0 D DualTec™ Kit (only in screw mode)
	CT LUCIA 602 <small>(Only available in the USA)</small>		+4.0 to +34.0 D R28 Injector / Z28 Cartridge

IOL-Injector Correlation Chart

Injector View



Injector	Cartridge		Injector Type	Incision Size
AT.Shooter A1-2000	VISCOJECT-BIO 1.8 Cartridge Set ^{1,2}		Reusable	1.8 mm
VISCOJECT-BIO 1.8 Injector Set ^{2,3}			Disposable	
Single Use A6 Injector Set			Disposable	
BLUEMIXS 180			Disposable, Preloaded	≥2.2 mm
AT.Shooter A1-2000	VISCOJECT-BIO 2.2 Cartridge Set ^{1,2}		Reusable	
VISCOJECT-BIO 2.2 Injector Set ²			Disposable	

IOL-Injector Correlation Chart

IOL View



Injector	Cartridge		Injector Type	Incision Size
BLUESERT 2.2 ⁵			Disposable	2.2 mm
BLUESERT 2.4			Disposable	2.4 mm
BLUESERT 2.6			Disposable	2.6 mm
SkyJet 2.8			Disposable, Preloaded	≥ 2.8 mm ≤ 3.0 mm
OD 655 DualTec kit 3-P Acrylic 2.8			Disposable	
R28 Injector ⁴	Z28 Cartridge ⁴		Reusable	
ACCUJECT 3.0			Disposable	

¹ Soft silicone cushion is included in the cartridge set

² VISCOJECT-BIO or "classic" versions

³ VISCOJECT 1.8 ≥ 1.8 mm with docking technique

⁴ Only available in the USA

⁵ Not sold as stand alone