

Sourcing Protective Coverings*

Materials, Dimensions and Potential Vendors*



Seeing beyond

	ZEISS Component	Protective Covering Type/Material	Protective Covering Dimensions		Vendor Name*	Vendor Item Number*
US	Forehead Rest Chinrest	Film	4 x 6 in	10 x 15 cm	Dental Savings Noble Dental Supplies Dental Office Products	BAR00-CLEAR 100-2102 D435-0002
	Keyboard	Bag, Polyethelene	7 x 18 in, 1 mil thick 4 x 8 in	20 x 40 cm, 50µm thick	U-Line International Plastics Inc	S-10877 PB00408
	Patient Response Button	Bag, Polyethelene	2 x 6 in, 2 mil thick 7 x 18 in	5 x 15 cm, 50µm thick	U-Line International Plastics Inc	S-1993 PB00718
	Forehead Rest	Bag	2 x 6 in, 2 mil thick 7 x 18 in	6.2 x 31 cm, 40µm thick	U-Line International Plastics Inc	S-3068 PB01014
	Mouse	Bag, Polyethelene	4 x 8 in, 1 mil thick 2 x 6 in	10 x 20, 25µm thick	U-Line International Plastics Inc	S-3450 PB20206
	Joystick 4.3 in/109 mm circumference	Bag, Polyethelene	3 x 4 in, 2 mil thick 3 x 4 in	6 x 8 cm, 50µm thick	U-Line International Plastics Inc	S-942 PB20304
	Keyboard 11.5 x 8.5 x 1 in thick (wired w/ glidepad below)	Bag, Polyethelene	10 x 14 in, 1 mil thick 3 x 14 in	25 x 35 cm, 25µm thick	U-Line International Plastics Inc	S-3203 PB20314
EU	Forehead Rest Chinrest	Film	10 x 15 cm		AliExpress	4000032722162
	Keyboard	Bag	25 x 40 cm 20 x 40 cm		Elke Plastic Grassholff	1FB250X40050 712840
	Patient Response Button	Bag	10 x 15 cm 10 x 15 cm		Elke Plastic Grassholff	1FB100X15050 712610
	Forehead Rest	Bag	6.2 x 31 cm		Elke Plastic	2FB062X31040
	Mouse	Bag	10 x 25 cm 12 x 180 cm		Elke Plastic Lacers Print & Packaging	1FB100X20025 1000-12018025
	Joystick 4.3 in/109 mm circumference	Bag	6 x 9 cm 6 x 8 cm		Elke Plastic Lacers Print & Packaging	1FB060X9050 1000-06008050
	Keyboard 11.5 x 8.5 x 1 in thick (wired w/ glidepad below)	Bag	25 x 35 cm 25 x 35 cm		Elke Plastic Lacers Print & Packaging	1FB250X35025 1000-25035025

*Film coverings are not sold by ZEISS. This information is provided to enable you to source protective coverings. Potential vendors are provided as an example only.

Sourcing Protective Coverings*

Materials, Dimensions and Potential Vendors*



Seeing beyond

	ZEISS Component	Protective Covering Type/Material	Protective Covering Dimensions	Vendor Name*	Vendor Item Number*
China	Forehead Rest Chinrest	Roll	10 cm x 15 cm	Kang Ming NA	N/A
	Keyboard	Plastic Bag	20 cm x 40 cm	Wang Nuo	N/A
	Patient Response Button	Plastic Bag	10 cm x 15 cm	Wang Nuo	N/A
	Forehead Rest	Plastic Bag	10 cm x 30 cm	Wang Nuo	N/A
	Mouse	Plastic Bag	14 cm x 20 cm	Wang Nuo	N/A
	Joystick 4.3 in/109 mm circumference	Plastic Bag	6 cm x 8 cm	Wang Nuo	N/A
	Keyboard 11.5 x 8.5 x 1 in thick (wired w/ glidepad below)	Plastic Bag	24 cm x 38 cm	Wang Nuo	N/A
India	Forehead Rest Chinrest	Film	W 4 x L 6	Neelkanth Healthcare Pvt. Ltd	N/A
	Keyboard	LDPE Bag	W 7 x L 18 x 25 microns	MICO Plast Industries Pvt. Ltd	N/A
	Patient Response Button	LDPE Bag	W 2 x L 6 x 50 microns	MICO Plast Industries Pvt. Ltd	N/A
	Forehead Rest	LDPE Bag	W 3 x L 14 x 50 microns	MICO Plast Industries Pvt. Ltd	N/A
	Mouse	LDPE Bag	W 4 x L 8 x 25 microns	MICO Plast Industries Pvt. Ltd	N/A
	Joystick 4.3 in/109 mm circumference	LDPE Bag	W 3 x L 4 x 50 microns	MICO Plast Industries Pvt. Ltd	N/A
	Keyboard 11.5 x 8.5 x 1 in thick (wired w/ glidepad below)	LDPE Bag	W 10 x L 14 x 25 microns	MICO Plast Industries Pvt. Ltd	N/A

*Film coverings are not sold by ZEISS. This information is provided to enable you to source protective coverings. Potential vendors are provided as an example only.
EN_38_200_00041 ©Carl Zeiss Meditec, Inc. 2020. All rights reserved.