


Detecting progression in the center of the visual field



Mary Durbin PhD¹, Thomas Callan OD¹, Sophia Yu¹, Robert Chang MD²,
Tin Aung PhD³, Ian P Conner MD, PhD⁴, Gary C Lee PhD¹

1. Carl Zeiss Meditec Inc, Dublin, CA, USA
2. Stanford University School of Medicine, Stanford, CA, USA
3. Singapore National Eye Centre, Singapore, Singapore, Singapore
4. University of Pittsburgh Medical Center, Pittsburgh, PA, USA



Disclosures

Carl Zeiss Meditec, Inc. (CZMI), Dublin, CA, USA

- Mary Durbin, PhD: CZMI (E)
- Gary Lee, PhD: CZMI (E)
- Thomas Callan, OD: CZMI (E)
- Sophia Yu, BS: CZMI (E)

Stanford University School of Medicine, Stanford, CA, USA

- Robert Chang, MD: CZMI (F)

Singapore National Eye Centre, Singapore, Singapore, Singapore

- Tin Aung, PhD: CZMI (F)


University of Pittsburgh Medical Center, Pittsburgh, PA, USA

- Ian P Conner, MD, PhD: CZMI (F); Ocugenix (S)





Background – HFA

Threshold Standard Automated Perimetry (SAP)



FORUM




*Differential Light Sensitivity



Purpose

- 24-2 visual field test has limited ability to detect change in the center of the visual field (VF).
- One approach to detecting change is to apply the change limits estimated from the central points in repeated 24-2 scans to 10-2 scans.
- In this study we use short-term reproducibility data to determine the specificity of using the 24-2 inner zone change limits from the HFA guided progression analysis (GPA) to detect change in the 10-2 VF.




Methods – Data for Test

- Three sites
- 74 eyes of 74 glaucoma subjects were enrolled
- HFA™ II-i (ZEISS, Dublin, CA)
- SITA Standard 10-2 and 24-2 VFs
- 5 repeat visits within 3 months
- Because the repeated visits occurred over a short period of time, we deemed any change observed in this data to be a false positive.




Methods – GPA

- In HFA, GPA change limits in the 24-2 field depend on:
 - baseline mean deviation (MD),
 - baseline pattern deviation (PD)
 - and location
- Two scans averaged to create baseline
- Follow ups compared to baseline
- Differences are compared to change limits by zones
- GPA is not performed where the MD or PD are too poor (censoring), denoted as "X"



Methods – Calculating False Positive Rate



- False positive rates (FPR) with 95% confidence intervals (CI) were pooled for all 74 patients and all 68 test points in the 10-2 for both follow-ups.
- Specificity defined as $1 - \text{FPR}$.
- FPR calculated two ways:
 - pooling all points,
 - pooling only uncensored points (those not marked “X”).



Methods / Results:

MD from 10-2 can be used in place of 24-2

- Correlation of 24-2 MD with 10-2 MD showed an R^2 of 0.67, a slope of 0.89 and an offset of 0.19, so it is reasonable to substitute the MD from 10-2 to determine change limits. (Consistent with Wu et al, *Ophthalmology Glaucoma* 2019)



Methods – change limits for 10-2 derived from 24-2 central zone

- Change limits as implemented in HFA Guided Progression Analysis
 - Limits = $f(\text{eccentricity, field status (MD), local defect depth (PD)})$


central zone


from 10-2 field


from 10-2 location

- Derived from 300+ eyes with 4 visits across multiple test strategies (but not 10-2)



Results – Specificity

Population Descriptive Statistics:

- Mean age was 63.6 (35.7 to 79.6) years
- Mean MD was -3.9 (-18.2 to 1.2) dB

False Positive Rate(s) / Specificity

- FPR was 0.047 (CI: 0.043 to 0.052) considering all 10064 pooled points
- FPR was 0.051 (CI: 0.046 to 0.055) considering 9380 pooled valid points
- Both FPRs consistent with specificity of ~95%. Similar performance for 24-2.
- Confirmation improves specificity to ~99%




Conclusions

- Specificity was ~95% for detecting progression in 10-2 fields using limits established in the central portion of the 24-2
- Simple GPA for 10-2 can be created using existing data
- Once limits have been established, sensitivity will depend on rate of change as compared to those limits



Detecting progression in the center of the visual field



Contact: mary.durbin@zeiss.com