

ZEISS CT SPHERIS 204

Technical Specifications



Seeing beyond



CT SPHERIS® 204

| | |
|---|--|
| Optic Design | Monofocal, spheric |
| Material | Hydrophilic acrylic (25 %) with hydrophobic surface properties |
| Optic Diameter | 6.0 mm |
| Total Diameter | 11.0 mm |
| Haptic Angulation | 0° |
| Lens Design | Single-piece |
| Incision Size | 2.2 mm |
| Company Labeled A-Constant ¹ | 118.2 |
| Diopter Range | From -10.0 to +30.0 D (from +31.0 to +45.0 D upon request) -10.0 to +10.0 D: 1.0 D increments +10.0 to +30.0 D: 0.5 D increments +30.0 to +45.0 D: 1.0 D increments |
| ACD | 5.01 |
| Implantation in | Capsular bag |
| Injector/Cartridge Set ² | For IOLs from -10.0 to +30.0 D*: AT.Shooter A1-2000/VISCOJECT™-BIO 2.2 Cartridge Set or VISCOJECT-BIO 2.2 Injector Set * For IOLs from +31.0 to +45.0 D ACCUJECT™ 3.0 Injector Set |

¹ For optimized A Constants and ACD Constants refer to IOLCon: www.iolcon.org

² Please refer to our web page for the most up-to-date references: www.zeiss.com/injectors

CE 0297

CT SPHERIS 204



Carl Zeiss Meditec AG

Goeschwitzer Strasse 51–52

07745 Jena

Germany

www.zeiss.com/iol

www.zeiss.com/med/contacts

EN_32_022_0097III Printed in Germany. CZ-IX/2019 International edition: Only for sale in selected countries

The contents of the brochure may differ from the current status of approval of the product or service offering in your country. Please contact our regional representatives for more information. Subject to changes in design and scope of delivery and due to ongoing technical development. CT SPHERIS is either a trademark or registered trademark of Carl Zeiss Meditec AG or other companies of the ZEISS Group in Germany and/or other countries. VISCOJECT is a trademark of Meditel (0482).
© Carl Zeiss Meditec AG, 2019. All rights reserved.