

ZEISS CT SPHERIS 209M

Technical Specifications



Seeing beyond



CT SPHERIS® 209M

Optic Design	Monofocal, spheric
Material	Hydrophilic acrylic (25 %) with hydrophobic surface properties
Optic Diameter	6.0 mm
Total Diameter	11.0 mm
Haptic Angulation	0°
Lens Design	Single-piece, MICS
Incision Size	1.8 mm
Company Labeled A-Constant ¹	117.9
Diopter Range	From 0.0 to +32.0 D 0.0 to +10.0 D, 1.0 D increments + 10.0 to +30.0 D, 0.5 D increments + 30.0 to +32.0 D, 1.0 D increments
ACD	4.82
Implantation in	Capsular Bag
Injector/Cartridge Set ²	AT.Shooter A1-2000/VISCOJECT™-BIO 1.8 Cartridge Set <i>or</i> VISCOJECT-BIO 1.8 Injector Set

¹ For optimized A Constants and ACD Constants refer to IOLCon: www.iolcon.org

² Please refer to our web page for the most up-to-date references: www.zeiss.com/injectors

CE 0297

CT SPHERIS 209M



Carl Zeiss Meditec AG

Goeschwitzer Strasse 51–52

07745 Jena

Germany

www.zeiss.com/iol

www.zeiss.com/med/contacts

EN_32_022_0098III Printed in Germany. CZ-VIII/2019 International edition: Only for sale in selected countries.
The contents of the brochure may differ from the current status of approval of the product or service offering in your country. Please contact our regional representatives for more information.
Subject to changes in design and scope of delivery and due to ongoing technical development. CT SPHERIS is either a trademark or registered trademark of Carl Zeiss Meditec AG or other companies of the ZEISS Group in Germany and/or other countries. VISCOJECT is a trademark of Mediceal (0482). A6 injector is a trademark of MDI (CE 0120).
© Carl Zeiss Meditec AG, 2019. All rights reserved.