

ZEISS COMBIVISC

Technical Specifications



Seeing beyond

	Cohesive part	Dispersive part
Behavior	Viscous-cohesive	Medium-viscosity dispersive
Origin	Bacterial fermentation	Bacterial fermentation
Substance	Sodium Hyaluronate	Sodium Hyaluronate
Concentration I (%)	1.5	3.0
Concentration II (mg/ml)	15	30
Volume (ml)	1.0	0.85
CDI	44	39
Osmolality (mOsmol/kg)	300 – 360	300 – 360
Molecular weight (megadalton ¹)	2.9	1.0
Pseudoplasticity	133	10
Zero-shear viscosity (Pa s ²)	295	47
Cannula	25 G	25 G
Storage temperature (°C)	2 – 8	2 – 8
Singlepack content:	2 x intracameral syringes, 2 x 25 G cannulas	
Multipack content:	10 x intracameral syringes, 10 x 25 G cannulas	

COMBIVISC® is latex and phthalate free, has a shelf life of 3 years and a pH range between 6.8 and 7.6.

¹ Megadalton = 1 million Dalton

² Pa s = 1000 mPa s



CE 0344
COMBIVISC



Carl Zeiss Meditec AG
Goeschwitzer Str. 51–52
07745 Jena
Germany
www.zeiss.com/ovd
www.zeiss.com/contacts



Hyaltech Ltd.
Starlaw Business Park
Livingston EH54 8SF
United Kingdom
www.zeiss.com/ovd
www.zeiss.com/contacts

EN_32_022_0221IV Printed in Germany. CZ-1/2021 International edition: Only for sale in selected countries.
The contents of the brochure may differ from the current status of approval of the product or service offering in your country. Please contact our regional representatives for more information. Subject to changes in design and scope of delivery and due to ongoing technical development.
COMBIVISC is a registered trademark of Carl Zeiss Meditec AG in Germany and/or other countries.
© Carl Zeiss Meditec AG, 2021. All rights reserved.