



Process No.	Title	Page
BLR-2-21-01-02-TP-08	IOL-Injector Correlation Chart	1 / 6

IOL		Diopter Range & Injector / Cartridge
4-Haptic MICS	CT SPHERIS® 209M	0.0 to +32.0 D AT.Shooter A1-2000 injector / VISCOJECT™-BIO 1.8 Cartridge or VISCOJECT-BIO 1.8 Injector Set or A6 injector & *AT.Smart cartridge
	CT ASPHINA® 409M	
	CT ASPHINA 509M	
	AT LISA® 809M	
	AT TORBI® 709M	-10.0 to +28.0 D sphere; +1.0 to +4.0 D cylinder AT.Shooter A1-2000 injector / VISCOJECT™-BIO 1.8 Cartridge or VISCOJECT-BIO 1.8 Injector Set or A6 injector & *AT.Smart cartridge -10.0 to +28.0 D sphere; +4.5 to +12.0 D cylinder VISCOJECT-BIO 2.2 Injector Set or AT.Shooter A1-2000 injector / VISCOJECT-BIO 2.2 cartridge Set +28.5 to +32.0 D sphere; +1.0 to +12.0 D cylinder VISCOJECT-BIO 2.2 Injector Set or AT.Shooter A1-2000 injector / VISCOJECT-BIO 2.2 cartridge Set
	AT LISA toric 909M	-10.0 to +24.0 D sphere; +1.0 to +12.0 D cylinder AT.Shooter A1-2000 injector / VISCOJECT™-BIO 1.8 Cartridge or VISCOJECT-BIO 1.8 Injector Set or A6 injector & *AT.Smart cartridge +24.5 to +32.0 D sphere; +1.0 to +12.0 D cylinder VISCOJECT-BIO 2.2 Injector Set or AT.Shooter A1-2000 injector / VISCOJECT-BIO 2.2 cartridge Set
	AT LISA tri toric 939M	-10.0 to +28.0 D sphere; +4.5 to +12.0 D cylinder and



Process No.	Title	Page
BLR-2-21-01-02-TP-08	IOL-Injector Correlation Chart	2 / 6

		<p>+28.5 to +32.0 D sphere; +1.0 to +12.0 D cylinder VISCOJECT-BIO 2.2 Injector Set or AT.Shooter A1-2000 injector / VISCOJECT-BIO 2.2 cartridge Set</p>
	AT LARA toric 929M	<p>-4.0 to +32.0 D SE; +4.5 to +12.0 D cylinder VISCOJECT-BIO 2.2 Injector Set or AT.Shooter A1-2000 injector / VISCOJECT-BIO 2.2 cartridge Set</p>
BLUEMIXS® Preloaded IOLs		
	CT ASPHINA 409MP	<p>0.0 to +32.0 D BLUEMIXS 180 injector</p>
	CT ASPHINA 509MP	
	AT LISA 809MP	
	AT LISA tri 839MP	
	AT LARA 829MP	
	AT TORBI 709MP	<p>-10.0 to +28.0 D sphere; +1.0 to +4.0 D cylinder BLUEMIXS 180 injector</p>
	AT LISA tri toric 939MP	<p>-10.0 to +28.0 D sphere; +1.0 to +4.0 D cylinder BLUEMIXS 180 injector</p>
	AT LARA toric 929MP	<p>-8.0D to +32.0 D SE; +1.0 to +4.0 D cylinder BLUEMIXS 180 injector</p>
4-Haptic	CT SPHERIS® 204	<p>-10.0 to +30.0 D AT.Shooter A1-2000 / VISCOJECT-BIO 2.2 Cartridge Set or VISCOJECT-BIO 2.2 Injector Set</p>
	CT ASPHINA® 404	<p>+31.0 to +45.0 D (CT SPHERIS 204) +31.0 to +42.0 D (CT ASPHINA 404) ACCUJECT™ 3.0 injector set</p>
3-Haptic	CT SPHERIS 203	<p>+8.0 to +32.0 D Forceps</p>



Process No.	Title	Page
BLR-2-21-01-02-TP-08	IOL-Injector Correlation Chart	3 / 6

<i>SkyJet® Preloaded IOLs</i>	
CT SPHERIS 203P	+8.0 to +26.0 D SkyJet 2.8 injector
	+27.0 to +30.0 D Forceps
CT ASPHINA 603P	+10.0 to +30.0 D SkyJet 2.8 injector
<i>Preloaded hydrophobic IOLs</i>	
CT LUCIA 601P/PY	+4.0 to +24.0 D ACCUJECT 2.0 injector set
	+24.5 to +30.0 D ACCUJECT 2.2 injector set
CT LUCIA 201P	+30.5 to +34.0 D ACCUJECT 2.4 injector set
	+4.0 to +24.0 D BLUEJECT 2.0 injector set
CT LUCIA 611P/PY	+24.5 to +30.0 D BLUEJECT 2.2 injector set
	+30.5 to +34.0 D BLUEJECT 2.4 injector set
CT LUCIA 211P	0.0 to +24.0 D BLUESERT 2.2 injector set
	+24.5 to +30.0 D BLUESERT 2.4 injector set
CT LUCIA 621P/PY	+30.5 to +34.0 D BLUESERT 2.6 injector set
	0.0 to +24.0 D
CT LUCIA 221P	0.0 to +24.0 D

1-Piece C-Loop



Process No.	Title	Page
BLR-2-21-01-02-TP-08	IOL-Injector Correlation Chart	4 / 6

		BLUESERT 2.2 injector set +24.5 to +30.0 D BLUESERT 2.4 injector set
3 Piece	CT LUCIA 202	+4.0 to +34.0 D OD655 DualTec™ Kit (only in screw mode)
	CT LUCIA 602 ⁴	+4.0 to +34.0 D R28 injector / Z28 cartridge



Process No.	Title	Page
BLR-2-21-01-02-TP-08	IOL-Injector Correlation Chart	5 / 6

Injector and Cartridge Type and Incision Size

Injector	Cartridge	Injector Type	Incision Size
AT.Shooter A1-2000	VISCOJECT-BIO 1.8 Cartridge set ^{1,2}	Reusable	1.8mm
VISCOJECT-BIO 1.8 Injector Set ²		Disposable	
Single Use Injector A6	*AT.Smart Cartridge Set ¹	Disposable	
BLUEMIXS 180		Disposable, Pre-Loaded	
AT.Shooter A1-2000	VISCOJECT-BIO 2.2 Cartridge set ^{1,2}	Reusable	≥2.2mm
VISCOJECT-BIO 2.2 Injector Set ²		Disposable	
ACCUJECT 2.0		Disposable	
ACCUJECT 2.2		Disposable	2.4mm
ACCUJECT 2.4		Disposable	2.6mm
BLUEJECT 2.0		Disposable	2.2mm
BLUEJECT 2.2		Disposable	2.4mm
BLUEJECT 2.4		Disposable	2.6mm
BLUESERT 2.2		Disposable	2.2mm
BLUESERT 2.4		Disposable	2.4mm
BLUESERT 2.6		Disposable	2.6mm
SkyJet 2.8		Disposable, Pre-Loaded	≥2.8 mm ≤ 3.0mm
OD 655 DualTec kit 3-P Acrylic 2.8		Disposable	
R28 injector ⁴	Z 28 cartridge ⁴	Reusable	
ACCUJECT 3.0		Disposable	



Process No.	Title	Page
BLR-2-21-01-02-TP-08	IOL-Injector Correlation Chart	6 / 6

--	--	--

¹ Soft silicone cushion is included in the cartridge set

² VISCOJECT - BIO or "classic" versions

³ VISCOJECT 1.8 ≥ 1.8mm with docking technique

⁴ Only available in the USA