

Safe Eyeglass Disinfection in 40 Seconds.

New
from
ZEISS



red dot winner 2021
industrial design



red dot winner 2021
innovative product



ZEISS UVClean.

The ZEISS UVClean disinfects **6 - 10 frames** in just **40 seconds** and kills up to **99.5%** of common viruses and bacteria on the surface of eyeglasses.*

www.zeiss.com/UVClean



Seeing beyond

Limited time offer
\$795
 Reg. Price \$995

Spend less time cleaning eyeglasses and more time selling them.

Don't settle for soapy water, smelly chemicals or re-purposed industrial/consumer UV solutions. The **ZEISS UVClean** was developed and rigorously tested to meet the demands of high volume eye care practices. The ZEISS UVClean is fast, safe and effective on frames & lenses.


1



Designed by ECPs to be stylish and fit in the dispensing area.


Easy to use with just **One Touch Start.**

3



ZEISS UVClean uses **4 germicidal UVC bulbs** to completely bathe each pair in powerful UV light.


2



Overall dimensions: 12" Wide x 12" High x 20" Deep and weighs 18 lbs.

Proprietary design specifically for eyeglasses accommodates **frames of all sizes.**

4



40 seconds to complete disinfection.
Kills up to 99.5% of common viruses and bacteria on the surface of eyeglasses.*

Why Choose ZEISS UVClean?

Designed specifically for ECPs ■ Proprietary 4 UVC Bulb System ■ 40 Second Cycle ■ No Liquids ■ No Harsh Chemicals ■ No Smells ■ No Disposable Waste ■ Small Footprint ■ Quiet Operation ■ One Touch Start ■ Safe & Gentle on Frames, Lenses & Frame Tags ■ ZEISS Branded Technology ■ The "Green Option" ■ **1-year warranty** ■ \$20 Shipping & Handling Fee ■ UL Listed ■ Made in the USA.

Limited number of pre-orders available on a first-come-first-served basis.

Reserve your **ZEISS UVClean** by contacting your **ZEISS** Business Development Representative or visit www.zeiss.com/UVClean today.

Carl Zeiss Vision Inc.

Customer Service USA: 1-866-596-5467

www.zeiss.com/UVClean

Follow us: Instagram: @ZEISSVisionCare_USA



Facebook: ZEISS Vision Care (US)



Seeing beyond

*Independent laboratory tests show effectiveness of 99.995% against adeno virus, 99.68% against influenza viruses, 99.47% against Escherichia coli, and 99.74% against staphylococcus aureus.
 ©2021 Carl Zeiss Vision Inc. 0000139.40417-ECP, Rev. 03/21