



## i.Terminal mobile

Start with digital centration – it's in your hands

Mobile centration is the new way to capture your patient's individual parameters with ZEISS high precision technology. Over 20 years of experience in centration systems are now available in its most convenient form yet, allowing for highest flexibility in the practice. The ZEISS i.Terminal mobile centration system for iPad gathers and calculates basic fitting parameters in a few seconds.

The application is designed to guarantee userfriendly control through intuitive handling to allow easy capturing of centration images. To further increase usability, the ZEISS i.Terminal mobile needs no bulky accessory to the iPad, such as additional cameras or flashes.

[Highlights](#)
[Technical Data](#)
[Details](#)
[Contact](#)

## Streamlined consultation is now combined in one solution




Your patient data is stored on a networking server of the portable i.Com mobile data management system by ZEISS in which the centration system is embedded. This allows for reliable procedures and easy data access from multiple iPads.

### Benefits:


- ZEISS i.Terminal mobile comes as part of ZEISS i.Com mobile, including the ZEISS i.Com server for storage of patient data and interfaces to common PMS.
- Simple and fast usage through face recognition.
- Smaller and lighter than a pupillometer.
- Including 2 calibration clips No additional equipment necessary.

# In 60 seconds, capture all necessary measuring parameters with i.Terminal mobile via iPad



- Frame data (A, B, DBL)
- Interpupillary distance (PD)
- Monocular pupillary distance (mono PD)
- Fitting height, segment height
- Back vertex distance (BVD)
- Pantoscopic angle (PA)
- Wrap angle


Interpupillary distance (PD)



Easy adjustment of fitting parameters by a fingertip on the iPad.




## Related Products



**i.Terminal 2**  
Another way to collect Centration Data

---

[Learn more](#)



**i.Com mobile**  
Streamline your workflow – capitalize your conversation flow

---

[Learn more](#)

ZEISS is one of the world's leading manufacturers of spectacle lenses, and is committed to delivering maximum precision and comfort. ZEISS designs and produces lenses, instruments and measurement systems, as well as retail concepts and technology services that continue to raise the bar in vision care.



**Contact us to get started**  
ZEISS Customer Service

[+49 7361 591-0](tel:+4973615910)

[Email](#)



Not all products, services or offers are approved or offered in every market and approved labelling and instructions may vary from one country to another. For country specific product information, see the appropriate country website.



---

**Explore**

- Understanding Vision
- Health + Prevention
- Lifestyle + Fashion
- Driving + Mobility
- Sports + Leisure
- Work Life

**Help me choose**

- Distance + Reading Glasses
- Varifocal Lenses
- Sunglasses
- Working Glasses
- Sports Glasses
- Glasses for Children
- Lens Coatings
- Contact Lenses
- Lens Cleaning
- At the Optician

**Services**

- My Vision Profile
- Online Vision Check

**For Eye Care Professionals**

- Keep Eyes Wide Open When Choosing Your Business Partner
- Instruments + Technologies
- ZEISS Spectacle Lenses
- ZEISS Cleaning Solutions

---

**We use cookies on this site.** Cookies are small text files that are stored on your computer by websites. Cookies are widely used and help to optimize the pages that you view. By using this site, you agree to their use. [> more](#)

OK