



## VISUREF® 100

The ZEISS VISUREF 100 is a 2-in-1 diagnostic instrument, including auto-refractometer and keratometer.

[Highlights](#)
[Details](#)
[Contact](#)


## 4 different operation modes

It provides accurate and reliable data for your daily practice routines. Intuitive handling features let you conduct the most important measurements with one device, including pupil and iris diameter evaluations. Refraction measurements of patients with implanted IOLs are also possible. Additionally, the automatic fogging incorporates an infinity scene to reduce the effects of accommodation.

### A valuable addition to your practice

When taking care of patients, every minute is valuable. There is no time for dealing with complex systems. Still, you need diagnostic data that is reliable and accurate, no matter what. Especially for the very first step in vision care: the refraction/keratometric evaluation.

**The ZEISS VISUREF 100. Adding more precision and reliability to your vision care – right from the start.**

#### Benefits:

- The motorised chinrest lets you quickly and comfortably position the patient



- Tilttable 6.4 inch TFT display to show diagnostic information
- Serial interface to export/transfer data to EMR/PMS or digital phoropter
- Quick release brake to disable XY movement of the instrument base and printout



### Tilttable 6.4 inch TFT display

Diagnostic information is clearly displayed in excellent contrast and resolution on the adjustable TFT monitor. Featuring a tilt range of 20°, it can be adjusted for individual user viewing preferences.



## Related Products



### Autorefraction with the i.Profiler plus

Fewer questions and more answers

[Learn more](#)

ZEISS is one of the world's leading manufacturers of spectacle lenses, and is committed to delivering maximum precision and comfort. ZEISS designs and produces lenses, instruments and measurement systems, as well as retail concepts and technology services that continue to raise the bar in vision care.

### Contact us to get started

ZEISS Customer Service

[+49 7361 591-0](#)

[Email](#)



Not all products, services or offers are approved or offered in every market and approved labelling and instructions may vary from one country to another. For country specific product information, see the appropriate country website.



---

**Explore**

Understanding Vision  
Health + Prevention  
Lifestyle + Fashion  
Driving + Mobility  
Sports + Leisure  
Work Life

**Help me choose**

Distance + Reading Glasses  
Varifocal Lenses  
Sunglasses  
Working Glasses  
Sports Glasses  
Glasses for Children  
Lens Coatings  
Contact Lenses

**Services**

My Vision Profile  
Online Vision Check  
Newsroom

**For Eye Care Professionals**

Our goal is to make you successful  
Instruments + Technologies  
ZEISS Spectacle Lenses  
ZEISS Cleaning Solutions

Lens Cleaning  
At the Optician  
Eye Surgery

---

[ZEISS International](#) | [Careers International](#) | [Press and Media International](#)

[Publisher](#) | [Legal Information](#) | [Data Protection](#) | [Cookie Preferences](#)