


## ZEISS Single Vision Superb

Widen your view

Beyond aspheric. Beyond standard atoric. The ZEISS Single Vision Superb lens for improved vision. Advanced freeform technology allows more precision compared to conventional single vision lenses, leading to clear and sharply focused vision over the entire lens.

[Design](#) [Benefit](#) [Details](#) [Product Availability](#)[Contact](#)

Point-by-point optimisation over the entire lens leads to improved vision

ZEISS Single Vision Superb lenses are designed for each individual wearer's prescription including: sphere, cylinder, axis, prism and prism base. Point-by-point optimisation is applied over the entire lens, leading to an even more comfortable single vision design.

# The benefits for the wearer

- ✓ Improved clarity and sharp vision
- ✓ Thinner & flatter lenses for improved cosmetic outcome
- ✓ Improved night and low light vision, if > [i.Scription®](#) technology is added

## ZEISS Single Vision Superb at a glance

- ✓ Using advanced freeform technology by ZEISS.
- ✓ Point-by-point optimisation of the lens surface leads to fewer power errors in the lens periphery and allows precise optical images.
- ✓ Optimised for every prescription: sphere, cylinder, axis, prism and prism base.
- ✓ Available with > [i.Scription®](#) technology.

## Product Availability

	Clear / tinted	PhotoFusion	Transitions®	Polarised
<b>Plastic 1.74</b>	■ (only clear)			
<b>Plastic 1.67</b>	■	■	■	
<b>Plastic 1.6</b> (also available in Sport)	■	■	■	■
<b>Polycarbonate 1.59</b> (also available in Sport)	■	■	■	■
<b>Trivex 1.53</b> (also available in Sport)	■	■	■	■
<b>Plastic 1.5</b> (also available in Sport)	■	■	■	■

ZEISS is one of the world's leading manufacturers of spectacle lenses, and is committed to delivering maximum precision and comfort. ZEISS designs and produces lenses, instruments and measurement systems, as well as retail concepts and technology

**Contact us to get started**  
ZEISS Customer Service

[+49 7361 591-0](#)

[Email](#)

services that continue to raise the bar in vision care.



Not all products, services or offers are approved or offered in every market and approved labelling and instructions may vary from one country to another. For country specific product information, see the appropriate country website.



---

**Explore**

Understanding Vision  
Health + Prevention  
Lifestyle + Fashion  
Driving + Mobility  
Sports + Leisure

**Help me choose**

Distance + Reading Glasses  
Varifocal Lenses  
Sunglasses  
Working Glasses  
Sports Glasses

**Services**

My Vision Profile  
Online Vision Check  
Newsroom

**For Eye Care Professionals**

Keep Eyes Wide Open When Choosing Your Business Partner  
Instruments + Technologies  
ZEISS Spectacle Lenses  
ZEISS Cleaning Solutions

Work Life

Glasses for Children

Lens Coatings

Contact Lenses

Lens Cleaning

At the Optician

---

[ZEISS International](#) | [Careers International](#) | [Press and Media International](#)

[Publisher](#) | [Legal Information](#) | [Data Protection](#) | [Cookie Preferences](#)