



ZEISS Sport Lens

Special lenses for fashionable wraparound frames and sunglasses.

ZEISS Sports Lenses provide the optimum optical correction for all curved frames - and makes aberrations a thing of the past.

ZEISS Sports Lenses - curved perfection

Whether for designer frames, highly fashionable sunglasses or sports eyewear – wrapped prescription lenses are definitely the latest trend. This confronts the optician with a special challenge: when frames of this type are glazed with traditional lenses, aberrations result for the user during wear.

The optimum solution


ZEISS Sports Lenses provide optimum optical correction for all wrapped frames - and it makes aberrations a thing of the past.

Benefits for the wearer

- ✓ Lightweight plastic lens for glazing highly curved frames (high bow angle)

- ✓ Individual optimization for optimum vision during wear
- ✓ High quality coating > [DuraVision® Platinum](#) available
- ✓ Available as single vision and as progressive lens
- ✓ Maximum diversity: same coatings possible as for > [ZEISS Single Vision Sph 1.5 / 1.6](#)
- ✓ ZEISS Single Vision Individual Sport is also available in > [photochromic](#) and > [polarized](#) versions

Wrapped Lenses with Cosmetic Edge® Technology



- ✓ ZEISS Cosmetic Edge Technology reduces the edge thickness in minus lenses by up to 30% in high prescription lenses
- ✓ Choice of frames: there's a fit for every frame thanks to very large lens diameters, variable decentration and flexible base curve selection.
- ✓ Excellent peripheral vision due to a fully-compensated optic design.
- ✓ Tailored lens designs from single vision to progressives, with great compensation for aberrations caused by lens tilt.
- ✓ Shaped to keep light and elements from entering the eye at all angles.

ZEISS Sport lenses with Cosmetic Edge® Technology Edge thickness reduced by up to 30%



Related Products



ZEISS Sun Sport Lenses

Wrapped lenses for sport frames offer clear optics and protection from wind and dust.

[Learn more](#)



ZEISS Lenses with UVProtect Technology

Sunglass-level UV protection in clear lenses.

[Learn more](#)



ZEISS DuraVision Platinum UV

The best of ZEISS Premium coatings

[Learn more](#)

manufacturers of eyeglass lenses, and is committed to delivering maximum precision and comfort. ZEISS designs and produces lenses, instruments and measurement systems, as well as retail concepts and technology services that continue to raise the bar in vision care.



ZEISS Customer Service

[1-866-596-5467](tel:1-866-596-5467)

[Email](#)

ZEISS Instruments Support

[1-800-201-4143](tel:1-800-201-4143)

Email: VTS-NA@zeiss.com



Explore

- [Understanding Vision](#)
- [Health + Prevention](#)
- [Lifestyle + Fashion](#)
- [Driving + Mobility](#)
- [Sports + Leisure](#)
- [Work Life](#)

Help me choose

- [Progressive Lenses](#)
- [Sunglasses](#)
- [Working Glasses](#)
- [Sports Glasses](#)
- [Glasses for Children](#)
- [Lens Coatings](#)
- [Lens Cleaning](#)
- [At the Eye Care Professional](#)

Services

- [My Vision Profile](#)
- [Online Vision Check](#)

For Eye Care Professionals

- [Our goal is to make you successful](#)
- [Instruments + Technologies](#)
- [ZEISS Eyeglass Lenses](#)
- [ZEISS Cleaning Solutions](#)